

**GI** *fx* **GI Effects**  
**Stool Profiles®**

# Interpretive Guide




# GI Effects Stool Profiles Interpretative Guide

The **GI Effects Interpretative Guide** has been created to provide a high-level approach to the GI Effects profile, biomarker interpretation, and therapeutic considerations. It is divided into two major sections: an overview of the GI Effects Interpretation At-a-Glance page; and a more in-depth review of the biomarkers comprising each of the Four Functional Pillars.

## Interpretation At-a-Glance Overview

Using evidence-based rules and weighted algorithms, the **Interpretation At-a-Glance** section on the first page of the GI Effects report synthesizes patient test results into key functional areas of clinical significance and provides a directional indication of potential next steps in patient management.



## Four Functional Pillars Biomarker Map

Infection Box	Inflammation Box	Insufficiency Box	Imbalance Box
<ul style="list-style-type: none"> <li>Any pathogenic organism present</li> </ul>	<ul style="list-style-type: none"> <li>Calprotectin</li> <li>Eosinophil Protein X (EPX)</li> <li>Fecal Secretory IgA</li> <li>Fecal Occult Blood</li> <li>Fecal Lactoferrin (if ordered)</li> </ul>	<ul style="list-style-type: none"> <li>Pancreatic Elastase 1 (PE1)</li> <li>Total Fecal Fats</li> <li>Products of Protein Breakdown (Total)</li> </ul>	<ul style="list-style-type: none"> <li>n-Butyrate</li> <li>Total SCFA</li> <li>Beta-glucuronidase</li> <li>Beneficial Bacteria <i>Lactobacillus</i>, <i>Bifidobacterium</i>, <i>E. coli</i> (PCR)</li> <li>Any potential pathogen (PP)</li> </ul>

## Four Functional Pillars

In this section, pertinent biomarkers have been grouped into four clinically actionable areas: Infection, Inflammation, Insufficiency, and Imbalance. The four functional pillars utilize a proprietary algorithm to evaluate key clinical markers in each of these four functional areas. The algorithm takes into account the level of each individual biomarker and its degree of clinical impact. As a result, an overall score of high, medium, or low is provided for each functional pillar. The score is represented by color-coded icons and informational graphics.

The specific biomarkers of concern that are utilized to establish the results for each functional pillar are listed above.

## Commensal Balance and Relative Abundance


The **Commensal Balance** infographic has been designed to provide a more precise view of an individual patient's commensal bacteria (PCR) results relative to a healthy cohort. It is a composite of two measures:

- The **Healthy-Pattern Continuum** (formerly known as the Diversity Association Index) is a progressive ranking scale based on a Genova proprietary algorithm that differentiates healthy and unhealthy commensal patterns. This algorithm is applied to an individual patient's GI Effects commensal bacteria (PCR) findings, and produces a numeric result ranging from 0 to 10 and is denoted by the 'y' axis of the Commensal Balance infographic.
- The **Reference Variance Score** reflects the total number of an individual patient's commensal bacteria (PCR) results that are out of reference range. This number ranges from zero to 24, and is denoted by the 'x' axis of the Commensal Balance infographic.

The patient's result on the Commensal Balance infographic is denoted by a black diamond against a color-coded gradient (green, yellow and red). The position of the patient's result against this background provides an At-a-Glance comparison of the patient's current commensal findings against those seen in healthy and diseased cohorts. Green suggests balanced commensal health status, yellow borderline, and red imbalanced.

The **Relative Abundance (RA)** graphic represents the proportional levels of selected phyla in an individual's microbiome and is represented relative to similar measures derived from a healthy cohort of individuals.


*See page 7 for commensal imbalance therapeutic recommendations.*



## FOUR FUNCTIONAL PILLARS BIOMARKER DETAIL

### INFECTION

This pillar is where common infectious microorganisms are reported and includes **pathogenic bacteria and intestinal parasites**.




### INFECTION



### INFLAMMATION


This pillar is where biomarkers that indicate inflammatory changes in the GI tract are reported. Biomarkers of GI inflammation and immunology provide information about the GI tract's interactions with, and responses to, the outside world. They indicate how well the GI tract is maintaining its role as a barrier, as well as whether the GI tract is undergoing pathological responses to external or internal challenges. The biomarkers are **Calprotectin**, a marker of neutrophil-driven inflammation<sup>2</sup>; **Eosinophil Protein X (EPX)**, a marker of eosinophil-driven inflammation and allergic response<sup>3</sup>; **Fecal Secretory IgA**, a marker of gut secretory immunity and barrier function<sup>4</sup>; and **Fecal Occult Blood**.




### INFLAMMATION




## EPX



## Fecal Secretory IgA



## Fecal Occult Blood




## INSUFFICIENCY

This pillar is where biomarkers that indicate digestive function are noted. Biomarkers of digestion and absorption provide information about nutrient breakdown and entry into the circulation. They ultimately indicate how well the GI tract is performing its basic digestive functions. The biomarkers are **Pancreatic Elastase 1**, a marker of exocrine pancreatic function<sup>12</sup>; **Products of Protein Breakdown**, markers of undigested protein reaching the colon<sup>13</sup>; and **Fecal Fat (Total)**, a marker of fat digestion/absorption.


## INSUFFICIENCY



### Pancreatic Elastase 1 (PE1)




### Products of Protein Breakdown (Total)



## Fecal Fat (Total)

---



### Target evaluation and treatment for common etiologies of fat maldigestion:

- **Pancreatic exocrine insufficiency**<sup>15</sup>
  - » If PE1 is less than 200, consider PERT
  - » Low fecal fat concentration does not exclude exocrine pancreatic insufficiency
- **Small Intestinal Bacterial Overgrowth (SIBO)**
  - » consider SIBO breath testing if: increased relative abundance, increased products of protein breakdown, increased SCFAs, or the presence of *Methanobrevibacter smithii*
- **Hypochlorhydria**<sup>16</sup>
  - » Assess for/reduce use of acid-blocking medications (as clinically indicated)
  - » Consider a betaine HCl challenge test, and treat as indicated
- **Bile Salt Insufficiency**<sup>15</sup>
  - » Assess for causes including liver damage, impaired gallbladder function
  - » Consider addition of bile salts and/or cholagogues

### Target evaluation and treatment for common conditions associated with fat malabsorption


- Infection<sup>17</sup>
- Celiac Disease<sup>18</sup>
  - » Consider Celiac and Gluten Sensitivity Panel
- IBS (confirm diagnosis via clinical criteria such as Rome IV)<sup>15,19</sup>
- IBD (review Calprotectin level; if greater than 120, GI referral)<sup>20</sup>
- Rapid transit time
- Gastric bypass, ileal resection or other surgeries that limit absorptive surface area<sup>21</sup>

### Further Evaluation:

- May be associated with deficiencies in fat or fat-soluble nutrients
  - » Consider nutritional assessment of essential fatty acids, fat-soluble vitamins



This pillar is where microbiome imbalances are noted. Biomarkers of the GI Microbiome provide information about the health, function, and abundance of the trillions of microbial cells in the GI tract. They indicate how well the microbiome is performing its shared metabolic functions with the human host. Abnormal results in this pillar may be associated with a large number of conditions and symptoms. The biomarkers which assess gut microbial imbalance are **Metabolic Indicators**, which demonstrates the microbiome's specific and vital metabolic function; these include **Short Chain Fatty Acids** that are produced by bacterial fermentation of fiber<sup>22</sup>, and **Beta-glucuronidase**, an enzyme produced by bacteria; **Commensal Bacteria**, a PCR evaluation of 24 key bacterial groups/species; and **Bacterial and Mycologic Culture**, which identifies potentially pathogenic (PP) organisms.<sup>23</sup>



## OVERALL PATTERN RECOGNITION

One of the benefits of ordering a comprehensive stool panel like GI Effects is the ability to observe trends across different groups of biomarkers. In addition to patterns discussed within the 4 functional pillars, clinicians commonly observe other patterns including Small Intestinal Bacterial Overgrowth (SIBO), Intestinal Permeability, Immune Dysregulation/ Loss of Resilience, and Inflammation. These patterns are not diagnostic, but rather suggestive of these common conditions and further workup may be appropriate.

Biomarker pattern, along with associated symptoms suggest	Biomarkers	Next steps
Small Intestinal Bacterial Overgrowth (SIBO)	<ul style="list-style-type: none"> <li>↑ Relative Abundance</li> <li>↑ Products of Protein Breakdown</li> <li>↑ SCFA</li> <li>↑ n-butyrate</li> <li>↑ Fecal Fat (total)</li> <li>↑ <i>Methanobrevibacter smithii</i></li> </ul>	<ul style="list-style-type: none"> <li>• Confirm with SIBO Breath Test</li> </ul>
Intestinal Permeability	<ul style="list-style-type: none"> <li>↑ sIgA</li> <li>↑ EPX</li> <li>↓ <i>Akkermansia muciniphila</i></li> </ul>	<ul style="list-style-type: none"> <li>• Confirm with Intestinal Permeability Assessment</li> <li>• Determine root cause of permeability</li> <li>• Consider GI repair/ support</li> <li>• Consider IgG and/or IgE Food Antibody testing</li> </ul>
Immune Dysregulation/ Loss of Resilience	<ul style="list-style-type: none"> <li>&lt;DL sIgA</li> <li>&lt;DL EPX</li> <li>↑ PP bacteria/yeast</li> <li>↓ Beneficial bacteria</li> <li>↓ <i>Bifidobacterium</i> spp.</li> <li>↑ <i>Methanobrevibacter smithii</i> and <i>Desulfovibrio piger</i></li> <li>High probability of <i>Blastocystis hominis</i></li> </ul>	<ul style="list-style-type: none"> <li>• Evaluate root cause of immune dysregulation (i.e., Adrenocortex Stress Profile)</li> <li>• Consider prebiotics/probiotics</li> <li>• Consider SIBO testing</li> </ul>

## References

1. Dinleyici EC, Eren M, Dogan N, Reyhanioglu S, Yargic ZA, Vandenplas Y. Clinical efficacy of *Saccharomyces boulardii* or metronidazole in symptomatic children with *Blastocystis hominis* infection. *Parasitology research*. 2011;108(3):541-545.
2. Steinbakk M, Naess-Andresen CF, Lingaas E, Dale I, Brandtzaeg P, Fagerhol MK. Antimicrobial actions of calcium binding leucocyte L1 protein, calprotectin. *Lancet*. 1990;336(8718):763-5.
3. Yang BG, Seoh JY, Jang MH. Regulatory Eosinophils in Inflammation and Metabolic Disorders. *Immune Netw*. 2017;17(1):41-7.
4. Mantis NJ, Rol N, Cortes B. Secretory IgA's complex roles in immunity and mucosal homeostasis in the gut. *Mucosal Immunol*. 2011;4(6):603-11.
5. D'Angelo F, Felley C, Frossard JL. Calprotectin in Daily Practice: Where Do We Stand in 2017? *Digestion*. 2017;95(4):293-301.
6. Majamaa H, Laine S, Miettinen A. Eosinophil protein X and eosinophil cationic protein as indicators of intestinal inflammation in infants with atopic eczema and food allergy. *Clin Exp Allergy*. 1999;29(11):1502-6.
7. van Odijk J, Peterson CG, Ahlstedt S, Bengtsson U, Borres MP, Hulthen L, et al. Measurements of eosinophil activation before and after food challenges in adults with food hypersensitivity. *Int Arch Allergy Immunol*. 2006;140(4):334-41.
8. Reimert CM, Tukahebwa EM, Kabatereine NB, Dunne DW, Vennervald BJ. Assessment of *Schistosoma mansoni* induced intestinal inflammation by means of eosinophil cationic protein, eosinophil protein X and myeloperoxidase before and after treatment with praziquantel. *Acta Trop*. 2008;105(3):253-9.
9. Bischoff SC, Grabowsky J, Manns MP. Quantification of inflammatory mediators in stool samples of patients with inflammatory bowel disorders and controls. *Dig Dis Sci*. 1997;42(2):394-403.
10. Wagner M, Sjoberg K, Vigren L, Olesen M, Benoni C, Toth E, et al. Elevated fecal levels of eosinophil granule proteins predict collagenous colitis in patients referred to colonoscopy due to chronic non-bloody diarrhea. *Scand J Gastroenterol*. 2016;51(7):835-41.
11. Brandtzaeg P. Update on mucosal immunoglobulin A in gastrointestinal disease. *Curr Opin Gastroenterol*. 2010;26(6):554-63.
12. Dominguez-Munoz JE, P DH, Lerch MM, Lohr MJ. Potential for Screening for Pancreatic Exocrine Insufficiency Using the Fecal Elastase-1 Test. *Digestive diseases and sciences*. 2017;62(5):1119-30.
13. Zarling EJ, Ruchim MA. Protein origin of the volatile fatty acids isobutyrate and isovalerate in human stool. *J Lab Clin Med*. 1987;109(5):566-70.
14. Ray A, Dittel BN. Interrelatedness between dysbiosis in the gut microbiota due to immunodeficiency and disease penetrance of colitis. *Immunology*. 2015;146(3):359-68.
15. Alkaade S, Vareedayah AA. A primer on exocrine pancreatic insufficiency, fat malabsorption, and fatty acid abnormalities. *Am J Manag Care*. 2017;23(12 Suppl):s203-s209.
16. Milovic V. Gastrointestinal disease and dosage form performance. In: *Oral Drug Absorption*. CRC Press; 2016:141-151.
17. Mohapatra S, Singh DP, Alcidi D, Pitchumoni CS. Beyond O&P Times Three. *The American journal of gastroenterology*. 2018;1.
18. Gujral N, Freeman HJ, Thomson ABR. Celiac disease: prevalence, diagnosis, pathogenesis and treatment. *World journal of gastroenterology*. 2012;18(42):6036-6059.
19. Watson L, Lalji A, Bodla S, Muls A, Andreyev HJ, Shaw C. Management of bile acid malabsorption using low-fat dietary interventions: a useful strategy applicable to some patients with diarrhoea-predominant irritable bowel syndrome? *Clinical medicine* (London, England). 2015;15(6):536-540.
20. Mourad FH, Barada KA, Saade NE. Impairment of Small Intestinal Function in Ulcerative Colitis: Role of Enteric Innervation. *Journal of Crohn's & colitis*. 2017;11(3):369-377.
21. O'Keefe SJD, Rakitt T, Ou J, et al. Pancreatic and Intestinal Function Post Roux-en-Y Gastric Bypass Surgery for Obesity. *Clinical and translational gastroenterology*. 2017;8(8):e112.
22. Rios-Covian D, Ruas-Madiedo P, Margolles A, Gueimonde M, de Los Reyes-Gavilan CG, Salazar N. Intestinal Short Chain Fatty Acids and their Link with Diet and Human Health. *Front Microbiol*. 2016;7:185.
23. Mroczynska M, Galecka M, Szachta P, Kamoda D, Libudzisz Z, Roszak D. Beta-glucuronidase and Beta-glucosidase activity in stool specimens of children with inflammatory bowel disease. *Polish journal of microbiology*. 2013;62(3):319-325.
24. Bach Knudsen KE. Microbial degradation of whole-grain complex carbohydrates and impact on short-chain fatty acids and health. *Adv Nutr*. 2015;6(2):206-213.
25. De Preter V, Raemen H, Cloetens L, Houben E, Rutgeerts P, Verbeke K. Effect of dietary intervention with different pre- and probiotics on intestinal bacterial enzyme activities. *Eur J Clin Nutr*. 2008;62(2):225-231.
26. Calcium-D-glucarate. *Altern Med Rev*. 2002;7(4):336-339.
27. Derrien M and Viegas P. Rethinking Diet to Aid Human-Microbe Symbiosis. *Trends in Microbiology*. 2017;25(2):100-112.
28. Le Chatelier E, et. al. Richness of human gut microbiome correlates with metabolic markers. *Nature*. 2013;500(7464):541-549.
29. Houghton D, Stewart C, Day C, Trenell M. Gut microbiota and lifestyle interventions in NAFLD. *Int J Mol Sci*. 2016;17(4):447.



Learn more at [www.gdx.net](http://www.gdx.net) or call 800.522.4762