



Oxidative Stress Analysis 2.0



Profile Components

RESERVE
Glutathione (GSH) Total Antioxidant Capacity (TAC) Cysteine (Cys-SH) Cystine (Cys-S-S-Cys) Cysteine/Cystine Ratio Sulphate Cysteine/Sulphate Ratio
PROTECTIVE ENZYMES
Superoxide Dismutase (SOD) Glutathione Peroxidase (GPx)
DAMAGE
Lipid Peroxides

What is Oxidative Stress?

Oxidative Stress occurs when the production of reactive oxygen species (ROS) outweighs the body's ability to remove them, thus shifting this equilibrium in the direction of oxidation (e.g. rusting). The instability of free radicals and other ROS causes them to extract electrons from neighboring molecules in a chain reaction, causing cellular damage in the process. Reducing agents, including dietary antioxidants, nutritional supplements, and antioxidant enzymes provide protection against free radical damage.

Oxidative stress has an integral relationship with the inflammatory cascade, which produces ROS, and is considered a driving force in the aging process. Oxidative stress has been implicated in a growing list of disorders, including cancer, atherosclerosis, arthritis, diabetes, macular degeneration, chronic fatigue syndrome, fibromyalgia, and neurodegenerative diseases.

Antioxidant needs vary significantly between individuals. Therefore, evaluating one's reduction/oxidation ("redox") balance can help pinpoint imbalances that are exacerbating or setting the stage for chronic illness, thus serving as a springboard for customized intervention. In this process we evaluate the individuals' unique biochemistry, including the availability of anti-oxidant RESERVES, the functionality of protective ENZYMES, and the presence (of absence) of tissue DAMAGE. Any imbalances are then treated specifically.

Overview of Pattern Interpretation:

Oxidative stress results from an imbalance between ROS production and antioxidant protection. Because Cystine is the oxidized form of Cysteine, a low Cysteine/Cystine Ratio can indicate a redox imbalance in the body, i.e., oxidative stress.


When prolonged or severe, oxidative stress eventually results in tissue damage and increased risk of disease, as indicated by an elevated Lipid Peroxides (reflecting oxidative damage to lipids in the body). The optional marker, 8-OHdG, can also shed light on oxidative damage to proteins.

The body's most important endogenous antioxidants include glutathione (GSH) and the enzymes, superoxide dismutase (SOD), glutathione peroxidase (GPx) and catalase. In the face of oxidative stress, the body generally compensates by increasing the synthesis of glutathione from component amino acids, and increasing the activity of antioxidant enzymes. Therefore, higher levels of GSH, SOD, or GPx suggests oxidative stress but also a healthy response to it. Normal Lipid Peroxides would suggest that despite the oxidative stress, oxidative damage has not yet occurred. However, it is still helpful to identify sources of oxidative stress and nutritionally protect these reserves.

Reserves—

Glutathione (GSH)

Glutathione (GSH) is a tripeptide made up of Glutamine, Glycine, and Cysteine and is the body's most potent endogenous anti-oxidant. GSH plays a central role in resistance to oxidative stress, functioning as an intracellular antioxidant, as well as a detoxifying agent for many xenobiotics. Low GSH contributes to oxidative stress and is thus a risk factor for many chronic diseases.




Low GSH associated with:

- Reduced antioxidant capacity, increased risk of oxidative damage and chronic illness
- Reduced ability to detoxify environmental toxins & metabolic byproducts
- Compromised gut lining, altered immunity, & reduced exercise endurance
- Reduced S-AdoMet synthesis & methylation

Higher GSH associated with:

- Efficient GSH synthesis
- Supplementation with GSH or its precursors
 - L-cysteine or
 - N-acetylcysteine (NAC),
 - Glycine,
 - L-glutamine)


Treatment Options to increase GSH



Reserves—

Total Antioxidant Capacity (TAC)


Antioxidants work synergistically, and wide individual variations exist regarding anti-oxidative capacity and needs. As a result, measuring total antioxidant response provides a more reliable indicator of antioxidant capacity than measuring single antioxidants. **Total Antioxidant Capacity (TAC)** reflects the collective power of reducing agents to neutralize free radicals for each individual.



Low TAC is associated with:

Oxidative stress and a need for greater antioxidant protection

Treatment Options to increase TAC



Cysteine

Cysteine is the rate-limiting amino acid for GSH, but also functions as an extracellular antioxidant and is a precursor for taurine, inorganic sulfate, acetyl-Coenzyme A, and protein synthesis.




Low Cysteine associated with:

- Reduced antioxidant capacity and ability to form GSH; increased risk of oxidative damage and chronic disease
- Possible low methionine (cysteine's precursor) or SAME (can be depleted during GSH upregulation)
- Possible high homocysteine (Hcy) and/or low B6 (converts Hcy to cysteine)
- Repeated DMSA administration (cysteine excretion) Possible cystinuria (inborn wasting of cysteine)

Higher Cysteine associated with:

- Supplementation with L-cysteine, N-acetylcysteine (NAC), L-methionine, or glutathione
- Possible blocks in cysteine metabolism (also see Cysteine/Sulfate Ratio, below)
- Possible cytotoxicity

Treatment Options to increase GSH



Reserves—

Cystine and 'Cysteine/Cystine Ratio'

Cystine is the oxidized disulfide form of cysteine (Cys) and is the predominant form of cysteine in the blood due to its greater relative stability. High cystine compared to cysteine, however, suggests a shifted redox balance and oxidative stress. The **Cys/CySS Ratio** is a reliable indicator of extracellular redox potential in the body. Low ratios have been noted in many chronic diseases.

Cystine (Cys-S-S-Cys)		1.60-3.22 mg/dL
Cysteine/Cystine Ratio		0.23-0.53


Low Cysteine/Cystine Ratio associated with:

- Shifted reduction/oxidation (redox) balance in the direction of oxidative stress
- Aging, smoking, increased risk of disorders such as type 2 diabetes, atherosclerosis, cancer

Higher Cysteine/Cystine Ratio associated with:

- Favorable 'redox' balance and antioxidant protection (see section for high Cysteine, if relevant)

Treatment Options to reduce Cystine and/or increase Cys/CySS Ratio



Sulfate and 'Cysteine/Sulfate Ratio'


Sulfate is produced from cysteine via sulfoxidation. Inorganic sulfate is a critical factor as part of Phase II detoxification reactions. Sulfate is also essential for mucin formation in the GI tract, and glycosaminoglycans in joint cartilage. Sulfoxidation defects result in low sulfate and increased risk of illness. The **Cysteine/Sulfate Ratio** reflects the efficiency of this 2-step conversion.

Sulfate		3.0-5.9 mg/dL
Cysteine/Sulfate Ratio		0.12-0.32

High Cysteine/Sulfate Ratio and/or low Sulfate is associated with:

- Leaky gut, environmental illness, food sensitivities, rheumatoid arthritis
- Alzheimer's, Parkinson's disease, motor neuron disease
- Sulfoxidation impairment (imbalance in enzymes, cysteine dioxygenase or sulfite oxidase)
- Low GSH (cysteine is used preferentially for GSH before being used for sulfate production)
 - Sulfite excess (can lower TAC and increase lipid peroxides)
 - Excessive exercise (can deplete GSH and sulfate)
 - Sulfate depletion from chronic arthritis or IBD
- Oxidative stress and a need for greater antioxidant protection

Treatment Options to increase TAC




Protective Enzymes—

Superoxide Dismutase & Glutathione Peroxidase

Superoxide Dismutase (SOD) and Glutathione Peroxidase (GPx) are endogenous antioxidant enzymes that protect against oxidative stress. SOD converts superoxide anion to hydrogen peroxide (H₂O₂), which is then inactivated by glutathione peroxidase (GPx) or catalase. SOD constitutes a critical step. However, if SOD is not balanced by GPx, H₂O₂ may still promote oxidative damage.


Superoxide Dismutase (SOD)	5,380	5,275-16,662 U/g Hb
Glutathione Peroxidase (GPx)	19.5	20.0-38.0 U/g Hb

Low SOD and Low GPx suggests:



Superoxide Dismutase (SOD)	15,967	5,275-16,662 U/g Hb
Glutathione Peroxidase (GPx)	36.4	20.0-38.0 U/g Hb


High SOD with high GPx suggests:



Protective Enzymes— Superoxide Dismutase & Glutathione Peroxidase (continued)


Superoxide Dismutase (SOD)		5,275-16,662 U/g Hb
Glutathione Peroxidase (GPX)		20.0-38.0 U/g Hb

High SOD with normal or low GPx suggests:



Superoxide Dismutase (SOD)		5,275-16,662 U/g Hb
Glutathione Peroxidase (GPX)		20.0-38.0 U/g Hb

High GPx with low/ normal SOD suggests:




Damage—

Lipid Peroxides

Lipid Peroxides are a direct indicator of oxidative damage to polyunsaturated fatty acids (PUFAs), suggesting that production of ROS has been inadequately balanced by antioxidants, esp. fat-soluble. Lipid peroxidation in cell membranes results in cellular dysfunction and is increased in many disease states. Lipid Peroxides are also increased in tissues that are poisoned by toxins.



High Lipid Peroxides are associated with:



This information is for the sole use of a licensed health care practitioner and is for educational purposes only. It is not meant for use as diagnostic information. All claims submitted to Medicare/Medicaid for Genova Diagnostics laboratory services must be for tests that are medically necessary. "Medically necessary" is defined as a test or procedure that is reasonable and necessary for the diagnosis or treatment of illness or injury or to improve the functioning of a malformed body member. Consequently, tests performed for screening purposes will not be reimbursed by the Medicare program.



Genova
Diagnostics®

Innovative Testing for Optimal Health

63 Zillicoa Street
Asheville, NC 28801
828.253.0621
www.GDX.net

© 2009 Genova Diagnostics
n,ig,oxstress2,110209