

Health Insights from Food and Environmental Antibody Testing



Antibody Assessments

CLINICIAN INFORMATION

GAIN INSIGHT INTO SYMPTOMS TRIGGERED BY ANTIGENIC SUBSTANCES

Genova Diagnostic's Food Antibody Assessments help identify those with true IgE-mediated allergies, as well as IgG-mediated food sensitivities. **The IgG Food Antibody Assessment** semi-quantitatively measures antibody levels to 87 foods and a total IgE level. **The IgE Food Antibody Profile** measures antibodies to the 19 most allergenic foods and a total IgE level. Additional tests are available for IgG vegetables, IgG spices, IgE regional inhalants, and IgE molds.

The key differences between IgE allergies and IgG sensitivities are summarized below:

IgE-Mediated Allergies (Foods, molds, inhalants)	IgG-Mediated Sensitivities (Foods, spices, vegetarian foods)
Immediate onset (minutes to hours)	Delayed onset (hours to days)
Circulating half-life of 1-2 days	Circulating half-life of 21 days
Permanent allergies	Temporary sensitivities
Stimulates histamine release	Activates complement Does not stimulate histamine release
Hives, stuffy or itchy nose, vomiting, diarrhea, wheezing, anaphylaxis, and other symptoms	Gastrointestinal symptoms, headaches, joint aches, rashes

Why Use Food Antibody Assessment?

Adverse food reactions can lead to distressing symptoms and chronic health conditions. Often times it is unknown exactly which food(s) may be the cause and testing can help identify the problematic foods. Removal of the reactive foods often results in resolution of symptoms.

Increased total antigenic load related to food and environmental reactivity has been associated with a wide range of medical conditions affecting virtually every part of the body. IgG-mediated sensitivities tend to be delayed and may be more vague or difficult to correlate with a specific food trigger.

Conditions associated with IgG food sensitivity:

- IBS¹⁻³
- Major Depressive Disorder³
- Migraine headaches⁴⁻⁶
- Skin rashes such as eczema⁷
- Joint aches⁸
- Autoimmune disease⁹
- Crohn's Disease¹⁰
- Obesity¹¹



**GASTROINTESTINAL
ISSUES**



MIGRAINES



ASTHMA

IgG Food Antibody Assessment (Serum)



63 Zillicoa Street
Asheville, NC 28801
© Genova Diagnostics

Patient: **SAMPLE PATIENT**

DOB:

Sex:

MRN:

IgG Food Antibody Results

Dairy	Vegetables	Fish/Shellfish	Nuts and Grains
Casein VL	Alfalfa 2+	Clam 0	Almond VL
Cheddar cheese 0	Asparagus 3+	Cod 1+	Buckwheat 3+
Cottage cheese 1+	Avocado 0	Crab 3+	Corn 1+
Cow's milk 1+	Beets 0	Lobster 3+	Corn gluten 2+
Goat's milk VL	Broccoli 1+	Oyster VL	Gluten VL
Lactalbumin 1+	Cabbage 1+	Red snapper 0	Kidney bean 1+
Yogurt VL	Carrot 0	Salmon 0	Lentil 2+
Fruits	Celery 0	Sardine 0	Lima bean 2+
Apple VL	Cucumber 1+	Shrimp 2+	Oat VL
Apricot 0	Garlic 0	Sole 0	Peanut VL
Banana 3+	Green Pepper VL	Trout 0	Pecan 2+
Blueberry 1+	Lettuce 1+	Tuna 0	Pinto bean VL
Cranberry 3+	Mushroom 1+	Poultry/Meats	Rice 1+
Grape VL	Olive 1+	Beef 0	Rye VL
Grapefruit 1+	Onion 0	Chicken 0	Sesame 3+
Lemon 2+	Pea VL	Egg white 1+	Soy 0
Orange 1+	Potato, sweet 0	Egg yolk 1+	Sunflower seed VL
Papaya 0	Potato, white 0	Lamb 0	Walnut 1+
Peach 0	Spinach VL	Pork 0	Wheat VL
Pear 0	String bean VL	Turkey 0	Miscellaneous
Pineapple 3+	Tomato VL		Yeast 1+
Plum VL	Zucchini VL		Cane sugar VL
Raspberry 0			Chocolate 1+
Strawberry 0			Coffee 1+
Total IgE			
Inside <input type="text"/> Outside <input type="text"/> Reference Range <input type="text"/>			
Total IgE ♦ <input type="text"/> 520.0 <=87.0 IU/mL			

0	None Detected	VL	Very Low	1+	Low	2+	Moderate	3+	High
---	---------------	----	----------	----	-----	----	----------	----	------

- The performance characteristics of all assays have been verified by Genova Diagnostics, Inc. Unless otherwise noted with ♦, the assay has not been cleared by the U.S. Food and Drug Administration.
- Total IgE level may have clinical significance regardless of specific antibody levels.
- Increasing levels of antibody detected suggest an increasing probability of clinical reactivity to specific foods.
- The Elimination Diet commentary is specific to IgG results only. Allergens inducing an IgE response should be completely avoided.

Laboratory Comments

Summary of IgG Test Results

Reactive / Non-Reactive Foods

3+ High

Asparagus	Banana	Buckwheat	Coconut
Crab	Cranberry	Curry	Garbanzo
Ginger	Lobster	Pineapple	Sesame
Vanilla			

2+ Moderate

Alfalfa	Bean sprout	Cashew	Corn gluten
Fennel	Lemon	Lentil	Lima bean
Oat bran	Pecan	Shrimp	Watermelon
Wild rice			

1+ Low

Blueberry	Broccoli	Cabbage	Chocolate
Cod	Coffee	Corn	Cottage cheese
Cow's milk	Cucumber	Cumin	Egg white
Egg yolk	Grapefruit	Kidney bean	Lactalbumin
Lettuce	Marjoram	Mushroom	Olive
Orange	Pistachio	Rice	Thyme
Walnut	Wheat bran	Yeast	

VL Very Low

Allspice	Almond	Apple	Basil
Black Pepper	Cane sugar	Cantaloupe	Casein
Cayenne	Cinnamon	Cloves	Filbert
Flax seed	Gluten	Goat's milk	Grape
Green pepper	Horseradish	Millet	Oat
Oyster	Paprika	Parmesan cheese	Pea
Peanut	Pinto bean	Plum	Rye
Sage	Spinach	String bean	Sunflower seed
Tomato	Wheat	Yogurt	Zucchini

0 None Detected

Apricot	Artichoke	Avocado	Bay leaf
Beef	Beets	Carrot	Celery
Cheddar cheese	Cherry	Chicken	Clam
Dill	Garlic	Kamut	Lamb
Mung bean	Mustard	Navy bean	Nutmeg
Onion	Oregano	Papaya	Parsley
Peach	Pear	Peppermint	Pork
Potato, sweet	Potato, white	Raspberry	Red Snapper
Rosemary	Safflower	Salmon	Sardine
Sole	Soy	Strawberry	Triticale
Trout	Tuna	Turkey	

Commentary

Overview

Immunoglobulin G (IgG) antibodies that elicit an immune response to food are in a class distinct from Immunoglobulin E (IgE) food allergy reactions. IgG-mediated food responses are described as delayed hypersensitivity reactions and have been associated in the peer-reviewed literature with an array of common clinical conditions including migraine, obesity, asthma, autoimmune diseases, and irritable bowel syndrome.

IgG Testing: Factors to Consider

IgG testing can be very useful in screening foods that a person is eating on a regular basis and which may be causing adverse reactions. However, it is possible to have adverse reactions to foods with low or non-detected levels of IgG. Because the IgG profile measures exposure of the immune system to food antigens, performing this test on a patient who is not consuming a particular food or who is taking a drug with known ability to suppress immune function (i.e. steroids) may result in the test not showing a positive reaction, potentially leading to a false negative result for the particular food. Be advised that if the patient is already on an elimination diet due to known food reactions, a negative result on an IgG food antibody profile does not necessarily mean that they can freely eat the food without experiencing symptoms.

IgG Results Interpretation

The amount of IgG antibodies is measured using a semi-quantitative ELISA assay procedure. The relative degrees of IgG present for each food are reported using a semi-quantitative level; None Detected (0), VL (very low), Low (1+), Moderate (2+) or High (3+). The degree of reactivity may not correlate with the severity of patient's response, therefore clinical correlation is advised as it can help direct treatment.

Clinical Management of Reactive IgG Foods: Elimination Diet

The purpose of an elimination diet is to pinpoint symptom-triggering foods that may be the root cause of and/or perpetuating chronic health issues. This diet is specific to food sensitivities that elicit an Immunoglobulin G (IgG) response and not those defined as classic (IgE-mediated) food allergy reactions. An elimination diet is a strategic process which depends on the oversight of the healthcare provider to ensure that a patient's nutritional requirements - macronutrient, micronutrient, and caloric needs - are adequate.

Four-Phases of an Elimination Diet



PHASE 1 – PREPARATION

A patient's clinical presentation and the IgG Food Antibody Assessment results typically determine which food(s) to temporarily remove from the diet. The average time frame for an elimination diet is 1 to 3 months. It is optimal to work with the patient to determine a start and end date for the elimination diet. Patient guidance around preparation ahead of the start date is important to ensure success. These include: (1) encouraging the patient to remove offending foods from the home and adjust grocery shopping accordingly; (2) providing the patient with resources that advance meal preparation, such as recipe books or reputable websites. Directing the patient to record foods consumed, date of consumption/elimination, and any notable changes in symptoms in a food journal can help track the progress of the diet.

Commentary



PHASE 2 – ELIMINATION

It is important to ensure the patient avoids those foods which resulted in a demonstrable reaction, either in whole food forms or as ingredients in other prepared foods to gain the greatest benefit. For patients unable to eliminate all reactive foods from their diet, focusing on the foods that elicited a stronger reaction (i.e.: 2+ and 3+) may be considered for an elimination diet. Practitioners may also encourage elimination of a complete food group when the patient shows reactivity to all foods tested within that group.



PHASE 3 – REINTRODUCTION

The reintroduction of eliminated foods is done one food at a time while monitoring for any adverse reaction. The patient should consume the test food several times throughout the day for several days. If symptoms occur with reintroduction, the patient should be instructed to remove that food once again and to evaluate whether the symptoms diminish over the next few days following elimination. Signs which may indicate an IgG food reaction include the following: headache, itching, bloating, fatigue, diarrhea or constipation, and indigestion. If the food does not cause symptoms during the reintroduction phase, it can be added back into the diet. The patient should continue this process with each food eliminated.

CAUTION: All patients warrant counseling related to signs and management of immediate hypersensitivity reactions prior to food reintroduction following an elimination diet. If reintroduction of a food causes an immediate allergic reaction (i.e. swelling of face, mouth, tongue, etc.; wheezing, rash/hives, or other allergic symptoms), it is imperative that the patient be treated as soon as possible. Following resolution of the immediate hypersensitivity reaction, the patient should be instructed to completely avoid consumption of that food.



PHASE 4 – LONG TERM MANAGEMENT

An elimination diet based on food sensitivity testing is part of a comprehensive approach to overall gastrointestinal health. Based on the test results and the complete clinical presentation of the patient, a long-term plan is usually developed utilizing the results of the reintroduction phase. Clinicians may also consider assessing and treating intestinal permeability, as gut barrier integrity is important for proper immune responses to foods. Nutrients that have been found to support intestinal barrier and decrease potential inflammation are glutamine, vitamin A, vitamin D, essential fatty acids (Omega-3), probiotics, and butyrate. Botanicals that can also be considered to assist with intestinal health are slippery elm, deglycyrrhizinated licorice (DGL), Aloe vera extract, and marshmallow root.

IgE Food Antibody Assessment



63 Zillicoa Street
Asheville, NC 28801
© Genova Diagnostics

Patient: **SAMPLE
PATIENT**

DOB:

Sex:

MRN:

IgE Food Antibody Results

	RESULT kU/L	CLASS	INDICATOR		RESULT kU/L	CLASS	INDICATOR
Grains				Nuts			
Buckwheat	<0.24	0/1	<input type="checkbox"/>	Almond	<0.24	0/1	<input type="checkbox"/>
Corn	<0.24	0/1	<input type="checkbox"/>	Brazil nut	<0.24	0/1	<input type="checkbox"/>
Oat	<0.24	0/1	<input type="checkbox"/>	Coconut	<0.24	0/1	<input type="checkbox"/>
Rice	<0.24	0/1	<input type="checkbox"/>	Hazelnut	<0.24	0/1	<input type="checkbox"/>
Sesame	<0.24	0/1	<input type="checkbox"/>	Peanut	<0.24	0/1	<input type="checkbox"/>
Soybean	<0.24	0/1	<input type="checkbox"/>	Seafood			
Wheat	<0.24	0/1	<input type="checkbox"/>	Blue mussel	<0.24	0/1	<input type="checkbox"/>
Dairy				Codfish	<0.24	0/1	<input type="checkbox"/>
Egg white	<0.24	0/1	<input type="checkbox"/>	Salmon	<0.24	0/1	<input type="checkbox"/>
Cow's milk	<0.24	0/1	<input type="checkbox"/>	Shrimp	<0.24	0/1	<input type="checkbox"/>
				Tuna	<0.24	0/1	<input type="checkbox"/>

Total IgE

	Inside	Outside	Reference Range
Total IgE	<input type="text" value=""/>	<input type="text" value="520.0"/>	<=87.0 IU/mL

- IgE levels must be used in conjunction with the clinical picture and are not intended to be independently diagnostic.
- The performance characteristics of all assays have been verified by Genova Diagnostics, Inc. All assays are cleared by the U.S. Food and Drug Administration.
- Total IgE level may have clinical significance regardless of specific antibody levels.
- Increasing levels of antibody detected suggest an increasing clinical reactivity to specific foods.

Key

Class	kU/L	Levels of Specific IgE	Indicator
0/1	<=0.24	Undetectable or Equivocal	<input type="checkbox"/>
I	0.25 - 0.39	Low	<input type="checkbox"/>
II	0.4 - 1.29	Moderate	<input type="checkbox"/>
III	1.3 - 3.89	High	<input type="checkbox"/>
IV	3.9 - 14.99	Very High	<input type="checkbox"/>
V	15 - 24.99	Very High	<input type="checkbox"/>
VI	>=25	Very High	<input type="checkbox"/>

Laboratory Comments

IgG Vegetarian Food Profile



63 Zillicoa Street
Asheville, NC 28801
© Genova Diagnostics

Patient: **SAMPLE PATIENT**

DOB:

Sex:

MRN:

<i>IgG Vegetable Food Results</i>					
Artichoke	0 <input type="checkbox"/>	Garbanzo	3+ <input type="checkbox"/>	Parmesan cheese	VL <input type="checkbox"/>
Bean sprout	2+ <input type="checkbox"/>	Filbert	VL <input type="checkbox"/>	Pistachio	1+ <input type="checkbox"/>
Cantaloupe	VL <input type="checkbox"/>	Kamut	0 <input type="checkbox"/>	Safflower	0 <input type="checkbox"/>
Cashew	2+ <input type="checkbox"/>	Millet	VL <input type="checkbox"/>	Triticale	0 <input type="checkbox"/>
Cherry	0 <input type="checkbox"/>	Mung bean	0 <input type="checkbox"/>	Watermelon	2+ <input type="checkbox"/>
Coconut	3+ <input type="checkbox"/>	Navy bean	0 <input type="checkbox"/>	Wheat bran	1+ <input type="checkbox"/>
Flax seed	VL <input type="checkbox"/>	Oat bran	2+ <input type="checkbox"/>	Wild rice	2+ <input type="checkbox"/>

<i>Total IgE</i>			
	Inside	Outside	Reference Range
Total IgE ♦	<input type="checkbox"/>	520.0	<=87.0 IU/mL

- The performance characteristics of all assays have been verified by Genova Diagnostics, Inc. Unless otherwise noted with ♦, the assays have not been cleared by the U.S. Food and Drug Administration.

- Increasing levels of antibody detected suggest an increasing probability of clinical reactivity to specific foods.

- Total IgE level may have clinical significance regardless of specific antibody levels.

0 <input type="checkbox"/>	None Detected	VL <input type="checkbox"/>	Very Low	1+ <input type="checkbox"/>	Low	2+ <input type="checkbox"/>	Moderate	3+ <input type="checkbox"/>	High
----------------------------	---------------	-----------------------------	----------	-----------------------------	-----	-----------------------------	----------	-----------------------------	------

Laboratory Comments



63 Zillicoa Street
Asheville, NC 28801
© Genova Diagnostics

Patient: **SAMPLE PATIENT**

DOB:

Sex:

MRN:

<i>IgG Spice Antibody Results</i>					
Allspice	VL		Curry	3+	
Basil	VL		Dill	0	
Bayleaf	0		Fennel	2+	
Black Pepper	VL		Ginger	3+	
Cayenne	VL		Horseradish	VL	
Cinnamon	VL		Marjoram	1+	
Cloves	VL		Mustard	0	
Cumin	1+		Nutmeg	0	
			Oregano	0	
			Paprika	VL	
			Parsley	0	
			Peppermint	0	
			Rosemary	0	
			Sage	VL	
			Thyme	1+	
			Vanilla	3+	

<i>Total IgE</i>			
	Inside	Outside	Reference Range
Total IgE ♦			<=87.0 IU/mL

- The performance characteristics of all assays have been verified by Genova Diagnostics, Inc. Unless otherwise noted with ♦, the assay has not been cleared by the U.S. Food and Drug Administration.

- Increasing levels of antibody detected suggest an increasing probability of clinical reactivity to specific foods.

- Total IgE level may have clinical significance regardless of specific antibody levels.

0		None Detected	VL		Very Low	1+		Low	2+		Moderate	3+		High
---	--	---------------	----	--	----------	----	--	-----	----	--	----------	----	--	------

Lab Comments

IgE Inhalants Profile

Texas +

63 Zillicoa Street
Asheville, NC 28801
© Genova Diagnostics




Patient: **SAMPLE PATIENT**

DOB:

Sex:

MRN:

IgE Antibody Levels			
INHALANT	RESULT	CLASS	INDICATOR
	kU/L		
Trees			
Maple	<0.24	0/1	<input type="checkbox"/>
Mountain Cedar	3.27	III	<input type="checkbox"/>
Grasses			
Bermuda Grass	0.67	II	<input type="checkbox"/>
June Grass (Kentucky Blue)	2.87	III	<input type="checkbox"/>
Perennial Rye Grass	3.57	III	<input type="checkbox"/>
Weeds			
Lamb's quarters	<0.24	0/1	<input type="checkbox"/>
English Plantain	<0.24	0/1	<input type="checkbox"/>
Rough Marsh Elder	<0.24	0/1	<input type="checkbox"/>
Giant Ragweed	<0.24	0/1	<input type="checkbox"/>
Molds			
Mold Generic	0.89	II	<input type="checkbox"/>
Misc.			
Cat dander	<0.24	0/1	<input type="checkbox"/>
Cockroach	<0.24	0/1	<input type="checkbox"/>
Dog dander	<0.24	0/1	<input type="checkbox"/>
Mite - D. farinae	0.57	II	<input type="checkbox"/>
Mite - D. microceras	0.77	II	<input type="checkbox"/>
Mite - D. pteronyssinus	0.41	II	<input type="checkbox"/>



- The performance characteristics of all assays have been verified by Genova Diagnostics, Inc. All assays are cleared by the U.S. Food and Drug Administration.
- Total IgE level may have clinical significance regardless of specific antibody levels.
- IgE levels must be used in conjunction with the clinical picture and are not intended to be independently diagnostic.

Total IgE		
	Inside	Outside
Total IgE	<input type="checkbox"/>	<input checked="" type="checkbox"/> 520.0
	Reference Range <=87.0 IU/mL	

Key			
Class	kU/L	Levels of Specific IgE	Indicator
0/1	<=0.24	Undetectable or Equivocal	<input type="checkbox"/>
I	0.25 - 0.39	Low	<input type="checkbox"/>
II	0.4 - 1.29	Moderate	<input type="checkbox"/>
III	1.3 - 3.89	High	<input type="checkbox"/>
IV	3.9 - 14.99	Very High	<input type="checkbox"/>
V	15 - 24.99	Very High	<input type="checkbox"/>
VI	>=25	Very High	<input type="checkbox"/>

Lab Comments

IgE Molds Profile



63 Zillicoa Street
Asheville, NC 28801
© Genova Diagnostics

Patient: **SAMPLE PATIENT**

DOB:

Sex:

MRN:

IgE Mold Antibody Results			
INHALANT	RESULT kU/L	CLASS	INDICATOR
Aspergillus fumigatus	<0.24	0/1	
Alternaria tenuis (Alternaria alternata)	3.12	III	
Candida albicans	<0.24	0/1	
Cladosporium herbarum	<0.24	0/1	
Curvularia lunata	0.36	I	
Epicoccum purpurascens	<0.24	0/1	
Fusarium moniliforme	<0.24	0/1	
Helminthosporium halodes	<0.24	0/1	
Mucor racemosus	<0.24	0/1	
Penicillium notatum	<0.24	0/1	
Phoma betae	0.4	II	
Pityrosporum orbiculare	0.42	II	
Rhizopus nigricans	0.53	II	
Stemphylium botryosum	0.81	II	
Trichoderma viride	0.25	I	

Key			
Class	kU/L	Levels of Specific IgE	Indicator
0/1	<=0.24	Undetectable or Equivocal	
I	0.25 - 0.39	Low	
II	0.4 - 1.29	Moderate	
III	1.3 - 3.89	High	
IV	3.9 - 14.99	Very High	
V	15 - 24.99	Very High	
VI	>=25	Very High	

- The performance characteristics of all assays have been verified by Genova Diagnostics, Inc. All assays are cleared by the U.S. Food and Drug Administration.
- Total IgE load may have clinical significance regardless of specific antibody levels.
- IgE levels must be used in conjunction with the clinical picture and are not intended to be independently diagnostic.

Total IgE		
	Inside	Outside
Total IgE		
		Reference Range <=87.0 IU/mL

Lab Comments



63 Zillicoa Street
Asheville, NC 28801
© Genova Diagnostics

Patient: **SAMPLE PATIENT**

DOB:

Sex:


MRN:

1006 Celiac Profile - Serum

Methodology: FEIA, Immunoturbidometric and IFA (when EMA IgA testing is performed)

Immunologic Markers			
Biomarker	Result		Reference Range
Total IgA	83	Insufficient	85-532 mg/dL
Anti-Tissue Transglutaminase IgG (tTG IgG)	2.7	Negative	<=6.9 U/ml
Anti-Deamidated Gliadin IgG (DGP IgG)	6.4	Negative	<=6.9 U/ml
Anti-Tissue Transglutaminase IgA (tTG IgA)	51.0	Positive	<=6.9 U/ml
Anti-Deamidated Gliadin IgA (DGP IgA)	6.4	Negative	<=6.9 U/ml
Anti-Endomysial IgA (EMA IgA)	Not Detected		Not Detected

Interpretation



● Related Profiles

Clinicians now have the ability, with a single requisition, to build any combination profiles using up to 7 antibodies profiles including:

See inside sample report for the full list of foods tested

- IgG Food Antibodies (87 IgG foods + total IgE) #1001
- IgG Spices (24 IgG spices + Total IgE) #1005
- IgG Vegetarian #1002
- IgE Food Antibodies (19 IgE Foods) #1000
- IgE Molds (15 IgE molds + Total IgE) #1004
- IgE Inhalants (14 IgE inhalants + Total IgE) #1003
- Celiac Profile #1018
- **A Full Combination of all seven (7) antibody profiles**

● References

1. Atkinson W, et. al. Food Elimination Based on IgG Antibodies in Irritable Bowel Syndrome: A Randomized Controlled Trial. *Gut*. 2004;53:1459-1464.
2. Drisko J, et. al. Treating Irritable Bowel Syndrome with a Food Elimination Diet Followed by Food Challenge and Probiotics. *J Amer Col Nutr*. 2006;25(6):514-522.
3. Karakula-Juchnowicz H, et.al. The food-specific serum IgG reactivity in major depressive disorder patients, irritable bowel syndrome patients and healthy controls. *Nutrients*. 2018 Apr;10(5):pii:E548.
4. Mitchell N, et. al. Randomized Controlled Trial of Food Elimination Diet Based on IgG Antibodies for the Prevention of Migraine-Like Headaches. *Nutr J*. 2011 May;10:85.
5. Alpay K, et. al. Diet Restriction in Migraine, Based on IgG Against Foods: A Clinical Double-Blind, Randomized, Cross-Over Trial. *Cephalgia*. 2010 Jul;30(7):829-37.
6. Aydinlar E, et.al. IgG-Based Elimination Diet in Migraine Plus Irritable Bowel Syndrome. *Headache*. 2013;53:514-525.
7. Shakib F, et. al. Study of IgG Sub-Class Antibodies in Patients with Milk Intolerance. *Clin Allergy*. 1986;16(5):451-458.
8. Panush RS. Food Induced ("Allergic") Arthritis: Clinical and Serologic Studies. *J Rheumatol*. 1990 Mar;17(3):291-4.
9. Coucke F. Food intolerance in patients with manifest autoimmunity. Observational study. *Autoimmune Rev*. 2018 Nov;17(11):1078-1080.
10. Bentz S, et. al. Clinical relevance of IgG antibodies against food antigens in Crohn's disease: a double-blind cross-over diet intervention study. *Digestion*. 2010;81(4):252-64.
11. Wilders-Truschning M, et. al. IgG antibodies against food antigens are correlated with inflammation and intima media thickness in obese juveniles. *Exp Clin Endocrinol Diabetes*. 2008 Apr;116(4):241-5.

