



**GI** *fx* **GI Effects**  
**Stool Profiles**



# The GI Effects® Advanced Interpretation

## Digging Deeper

Michelle Maddux, ND  
Senior Manager, Medical Education  
Genova Diagnostics

**GENOVA**  
DIAGNOSTICS®



# Lahnor Powell, ND, MPH

Medical Education Specialist - Atlanta




**Michelle Maddux, ND**  
Senior Manager, Medical Education



# Technical Issues & Clinical Questions

- Please type any technical issue or clinical question into either the “Chat” or “Questions” boxes, making sure to send them to “Organizer” at any time during the webinar.
- We will be compiling your clinical questions and answering as many as we can the final 15 minutes of the webinar.



DISCLAIMER: Please note that any and all emails provided may be used for follow up correspondence and/or for further communication.

**Need More Resources?**  
*Ensure you have an account!*

The screenshot shows the Genova Diagnostics website homepage. At the top is a navigation bar with social media icons and links for International, About Us, Contact Us, Search, myGDX, and US. Below this is the Genova Diagnostics logo and a main navigation menu with Home, Clinicians, and Patients. A banner image of a man and woman in a kitchen is featured, with the text "NutrEval® with Genomics" and the tagline "The Nutritional Test You Rely On Just Got Better!". A "LEARN MORE" button is present. Below the banner are three white boxes on a blue background, each with a red circle around it: "Getting Started" (with a hand icon), "Test Menu" (with a magnifying glass icon), and "MYGDX Login" (with a lock icon). At the bottom, another red circle highlights the "Online Education" section, which includes a description of free resources and a "LEARN NOW" button.

GENOVA  
DIAGNOSTICS

International About Us Contact Us Search myGDX US

HOME CLINICIANS PATIENTS

Website Preview - PLEASE SHARE YOUR FEEDBACK

**NutrEval® with Genomics**  
*The Nutritional Test You Rely On Just Got Better!*

LEARN MORE »

**Getting Started**  
Simple account setup. Licensed healthcare practitioners may begin the process of opening a free account here.  
NEW USERS

**Test Menu**  
A comprehensive menu of our diagnostic tests, including test descriptions, specimen requirements and kit instructions.  
SEARCH TESTS

**MYGDX Login**  
Clinicians: Log in to the MYGDX™ portal to order test kits and materials, download patient results, edit account information.  
LOG IN

**Online Education**  
Visit our Medical Education section for access to LIVEGDX Webinars, Educational Modules, Conferences, and LEARN GDX – short learning modules that demonstrate the clinical utility and diagnostic significance of key biomarkers. The modules are absolutely free to view!  
LEARN NOW

Start Using These Free Resources Today



**GI** *fx* **GI Effects**  
**Stool Profiles**



# The GI Effects® Advanced Interpretation

Digging Deeper

Michelle Maddux, ND  
Senior Manager, Medical Education  
Genova Diagnostics

**GENOVA**  
DIAGNOSTICS®

# Learning Objectives

- Gain understanding of the complexity of the microbiome and its metabolic importance
- Recognize physiologic digestive and absorptive patterns and inflammation seen in specific clinical conditions
- Apply the GI Effects® Comprehensive Stool Profile in clinical practice






# Case Example

- 52 y/o female with bloating, excessive belching and flatulence
- Constipation
  - Uses a laxative three times per week to achieve a bowel movement
- Cerebellar ataxia
- Overweight
- Sedentary lifestyle
- WNL thyroid function
- Eats a standard American diet; sugar cravings & avoids fiber (aggravates symptoms)



# Interpretation-at-a-Glance: Developed using Evidence-based Rules and Weighted Algorithms


- Functional Pillars of Health
  - Infection
  - Inflammation
  - Insufficiency (Digestive)
  - Imbalance (Metabolic)
- Global Gut Health Markers
  - Commensal Balance
  - Relative Abundance



# The Four Functional Pillars:


Provide immediate insight into clinically actionable findings in the areas of:

- Infection
- Inflammation
- Insufficiency
- Imbalance






# Commensal Balance and Relative Abundance Graphics Serve as General Markers of Gut Health



**Commensal Balance (CB):** Healthy-Pattern Continuum & Reference Variance Score  
**Relative Abundance (RA):** Identifies the level of growth of select bacterial phyla



- The **Commensal Balance** is a composite of two measures:
  - ‘Y-axis’: The **Healthy-Pattern Continuum** (formerly known as the Diversity Association Index) is a progressive ranking scale which differentiates healthy and unhealthy commensal patterns
  - ‘X-axis’: The **Reference Variance Score** reflects the total number of an individual patient’s commensal bacteria (PCR) results that are out of reference range





# Commensal Balance


- Intersection of two variables – The Healthy Pattern Continuum and the Reference Variance Score - is denoted by black diamond against a color-coded gradient (green, yellow and red)
- The position of the patient's result against this background provides an *At-a-Glance* comparison of the patient's current commensal findings against those seen in healthy and diseased cohorts
  - Green suggests balanced commensal health status
    - Represents 95% of individuals in the healthy cohort
  - Yellow suggests borderline commensal health status
  - Red suggests imbalanced commensal health status



# Relative Abundance

Identifies the level of growth of each measured bacterial phyla

- Significantly lower total Relative Abundance may represent an imbalance due to low beneficial commensal organisms
- Significantly higher total RA may represent an imbalance due to microbial overgrowth



# Commensal and Biomarker Clinical Association Charts

Patient: JANE DOE

ID:

Page 2

2200 GI Effects™ Comprehensive Profile – Stool

Interpretation At-a-Glance

Commensal Bacteria	Patient Results Out of Reference Range	Genova Diagnostics Commensal Bacteria Clinical Associations*							
		IBS	IBD	Metabolic Syndrome	Chronic Fatigue	Auto-immune	Type 2 Diabetes	High Blood Pressure	Mood Disorders
<b>Bacteroidetes Phylum</b>									
<i>Bacteroides-Prevotella</i> group		↑	↑	↑	↑	↑	↑	↑	↑
<i>Bacteroides vulgatus</i>	H	↑			↑	↑		↑	↑
<i>Barnesiella</i> spp.	H								
<i>Odoribacter</i> spp.	H								
<i>Prevotella</i> spp.	H	↑		↑	↑	↑		↑	↑
<b>Firmicutes Phylum</b>									
<i>Anaerotruncus colihominis</i>	H	↑	↑	↑	↑	↑	↑	↑	↑
<i>Butyrivibrio crossotus</i>									
<i>Clostridium</i> spp.	H								
<i>Coproccoccus eufactus</i>	H	↑			↑	↑		↑	↑
<i>Faecalibacterium prausnitzii</i>	H	↑				↑			↑
<i>Lactobacillus</i> spp.	H								
<i>Pseudoflavonifractor</i> spp.	H	↑	↑	↑	↑	↑	↑	↑	↑
<i>Roseburia</i> spp.			↓						
<i>Ruminococcus</i> spp.	L	↑↓	↓	↓	↓	↑↓	↑↓	↑↓	↑↓
<i>Veillonella</i> spp.	H	↑	↑	↑	↑	↑	↑	↑	↑
<b>Actinobacteria Phylum</b>									
<i>Bifidobacterium</i> spp.									
<i>Bifidobacterium longum</i>									
<i>Collinsella aerofaciens</i>		↑↑	↑↑	↓	↑↑	↑↑	↑↑	↑↑	↑↑
<b>Proteobacteria Phylum</b>									
<i>Desulfovibrio piger</i>									↑
<i>Escherichia coli</i>	H	↑	↑	↑	↑	↑	↑	↑	↑
<i>Oxalobacter formigenes</i>	H	↑		↑	↑				↑
<b>Euryarchaeota Phylum</b>									
<i>Methanobrevibacter smithii</i>	H	↑				↑			↑
<b>Fusobacteria Phylum</b>									
<i>Fusobacterium</i> spp.		↑	↑	↑	↑	↑	↑	↑	↑
<b>Verrucomicrobia Phylum</b>									
<i>Akkermansia muciniphila</i>		↓	↓	↓	↓	↓	↓	↓	↓
*Information derived from GDx results data comparing a healthy cohort to various clinical condition cohorts. The chart above showing a comparison of patient results to clinical conditions is meant for informational purposes only; it is not diagnostic, nor does it imply that the patient has a specific clinical diagnosis or condition.									
The arrows indicate Genova's clinical condition cohort test results falling below ↓ or above ↑ the reference range that is greater than that of Genova's healthy cohort.									
↑↓ Indicates Genova's clinical condition cohort test results falling below and above the reference range that are greater than that of Genova's healthy cohort.									
Cells with bolded arrows indicate Genova's clinical condition cohort had more test results falling above versus below ↑↓ or more below versus above ↓↑ the reference range compared to that of Genova's healthy cohort.									

© Genova Diagnostics - Robert M. David, PhD, Lab Director - CLIA Lic. #11D0255349 - Medicare Lic. #34-8475 - Georgia Lab Lic. Code #067-007  
 New York Clinical Lab PFI #4578 - Florida Clinical Lab Lic. #800008124

Patient: JANE DOE

ID:

Page 3

2200 GI Effects™ Comprehensive Profile – Stool

Interpretation At-a-Glance

Biomarker	Patient Results Out of Reference Range	Genova Diagnostics Biomarker Clinical Associations*							
		IBS	IBD	Metabolic Syndrome	Chronic Fatigue	Auto-immune	Type 2 Diabetes	High Blood Pressure	Mood Disorders
Pancreatic Elastase		↓	↓	↓	↓	↓	↓	↓	↓
Products of Protein Breakdown (Total)	L						↓↑		
Fecal Fat (Total*)	H	↑		↑	↑	↑	↓↑	↑	↑
Triglycerides		↑			↑	↑	↑	↑	↑
Long Chain Fatty Acids	H	↑			↑	↑	↓↑	↑	↑
Cholesterol							↓↑	↑	
Phospholipids	H	↑	↑	↑	↑	↑	↑	↑	↑
Calprotectin			↑					↑	
Eosinophil Protein X (EPX)			↑						
Fecal sIgA		↑	↑	↑	↑	↑	↑	↑	↑
Short Chain Fatty Acids (SCFA) (Total)	L				↓	↓			
n-Butyrate Concentration	L			↓					
n-Butyrate %	H								
Acetate%	L				↓↑		↓↑		
Propionate %				↑			↑	↑	
Beta-glucuronidase						↓↑			↓↑

\*Information derived from GDx results data comparing a healthy cohort to various clinical condition cohorts. The chart above showing a comparison of patient results to clinical conditions is meant for informational purposes only; it is not diagnostic, nor does it imply that the patient has a specific clinical diagnosis or condition.

The arrows indicate Genova's clinical condition cohort test results falling below ↓ or above ↑ the reference range that is greater than that of Genova's healthy cohort.

↓↑ Indicates Genova's clinical condition cohort test results falling below and above the reference range that are greater than that of Genova's healthy cohort.

Cells with bolded arrows indicate Genova's clinical condition cohort had more test results falling above versus below ↓↑ or more below versus above ↑↓ the reference range compared to that of Genova's healthy cohort.

© Genova Diagnostics - Robert M. David, PhD, Lab Director - CLIA Lic. #11D0255349 - Medicare Lic. #34-8475 - Georgia Lab Lic. Code #067-007  
New York Clinical Lab PFI #4578 - Florida Clinical Lab Lic. #600008124



# Commensal and Biomarker Clinical Association Charts

- Based on Genova Diagnostics GI Effects test-results database; the Commensal and Biomarker Clinical Association charts were developed
  - Based on self-reported clinical conditions (IBS, IBD, Metabolic Syndrome, Chronic Fatigue, Autoimmune dysfunction, Type 2 Diabetes, High Blood Pressure, and Mood Disorders)
- Differences between the healthy cohort and individuals with clinical conditions are denoted by the arrows in the Clinical Association charts

The arrows indicate Genova's clinical condition cohort test results falling below ↓ or above ↑ the reference range that is greater than that of Genova's healthy cohort.

↑↓ Indicates Genova's clinical condition cohort test results falling below and above the reference range that are greater than that of Genova's healthy cohort.

Cells with bolded arrows indicate Genova's clinical condition cohort had more test results falling above versus below ↑↓ or more below versus above ↓↑ the reference range compared to that of Genova's healthy cohort.



# Commensal Bacteria Clinical Association Chart

Patient: JANE DOE ID: Page 2

2200 GI Effects™ Comprehensive Profile – Stool

Interpretation At-a-Glance

Commensal Bacteria	Patient Results Out of Reference Range	IBS	IBD	Metabolic Syndrome	Chronic Fatigue	Auto-immune	Type 2 Diabetes	High Blood Pressure	Mood Disorders
<b>Bacteroidetes Phylum</b>									
Bacteroides-Prevotella group		↑	↑	↑	↑	↑	↑	↑	↑
Bacteroides vulgatus	H	↑			↑	↑		↑	↑
Barnesiella spp.	H								
Oribacter spp.	H								
Prevotella spp.	H	↑		↑	↑	↑		↑	↑
<b>Firmicutes Phylum</b>									
Anaerotruncus colihominis	H	↑	↑	↑	↑	↑	↑	↑	↑
Butyrivibrio crossotus									
Clostridium spp.	H								
Coprococcus eutactus	H	↑			↑	↑		↑	↑
Faecalibacterium prausnitzii	H	↑				↑			↑
Lactobacillus spp.	H								
Pseudoflavonifractor spp.	H	↑	↑	↑	↑	↑	↑	↑	↑
Roseburia spp.			↓						
Ruminococcus spp.	L	↑	↓	↓	↓	↑	↑	↑	↑
Veillonella spp.	H	↑	↑	↑	↑	↑	↑	↑	↑
<b>Actinobacteria Phylum</b>									
Bifidobacterium spp.									
Bifidobacterium longum									
Collinsella aerofaciens		↑	↑	↓	↑	↑	↑	↑	↑
<b>Proteobacteria Phylum</b>									
Desulfovibrio piger									↑
Escherichia coli	H	↑	↑	↑	↑	↑	↑	↑	↑
Oxalobacter formigenes	H	↑		↑	↑				↑
<b>Euryarchaeota Phylum</b>									
Methanobrevibacter smithii	H	↑				↑			↑
<b>Fusobacteria Phylum</b>									
Fusobacterium spp.		↑	↑	↑	↑	↑	↑	↑	↑
<b>Verrucomicrobia Phylum</b>									
Akkermansia muciniphila		↓	↓	↓	↓	↓	↓	↓	↓

\*Information derived from GDX results data comparing a healthy cohort to various clinical condition cohorts. The chart above showing a comparison of patient results to clinical conditions is meant for informational purposes only, it is not diagnostic, nor does it imply that the patient has a specific clinical diagnosis or condition.

The arrows indicate Genova's clinical condition cohort test results falling below ↓ or above ↑ the reference range that is greater than that of Genova's healthy cohort.

↑ ↓ Indicates Genova's clinical condition cohort test results falling below and above the reference range that are greater than that of Genova's healthy cohort.

Cells with bolded arrows indicate Genova's clinical condition cohort had more test results falling above below ↓ or more below versus above ↑ the reference range compared to that of Genova's healthy cohort.

© Genova Diagnostics · Robert M. David, PhD, Lab Director · CLIA Lic. #11D0255349 · Medicare Lic. #34-8475 · Georgia Lab Lic. Code #067-007  
New York Clinical Lab PFI #4576 · Florida Clinical Lab Lic. #60008124

- Detailed Commensal Bacteria (PCR) patient results on a grid that also contains the details for GDX microbial clinical associations for various health conditions
  - Allows direct comparison of the patient's current commensal patterns to aggregate GDX commensal patterns seen in both healthy and unhealthy
  - “apples-to-apples” comparison provides more relevant information for results interpretation of commensal-bacteria status



# Biomarker Clinical Association Chart

Patient: JANE DOE ID: Page 3

2200 GI Effects™ Comprehensive Profile – Stool

Biomarker	Patient Results Out of Reference Range	Genova Diagnostics Biomarker Clinical Associations*							
		IBS	IBD	Metabolic Syndrome	Chronic Fatigue	Auto-immune	Type 2 Diabetes	High Blood Pressure	Mood Disorders
Pancreatic Elastase		↓	↓	↓	↓	↓	↓	↓	↓
Products of Protein Breakdown (Total)	L						↓↑		
Fecal Fat (Total*)	H	↑		↑	↑	↑	↓↑	↑	↑
Triglycerides		↑			↑	↑	↑	↑	↑
Long Chain Fatty Acids	H	↑			↑	↑	↓↑	↑	↑
Cholesterol							↓↑	↑	
Phospholipids	H	↑	↑	↑	↑	↑	↑	↑	↑
Calprotectin			↑					↑	
Eosinophil Protein X (EPX)			↑						
Fecal sIgA		↑	↑	↑	↑	↑	↑	↑	↑
Short Chain Fatty Acids (SCFA) (Total)	L				↓	↓			
n-Butyrate Concentration	L			↓					
n-Butyrate %	H								
Acetate%	L				↓↑		↓↑		
Propionate %				↑			↑	↑	
Beta-glucuronidase						↓↑			↓↑

\*Information derived from GDx results data comparing a healthy cohort to various clinical condition cohorts. The chart above showing a comparison of patient results to clinical conditions is meant for informational purposes only; it is not diagnostic, nor does it imply that the patient has a specific clinical diagnosis or condition.

The arrows indicate Genova's clinical condition cohort test results falling below ↓ or above ↑ the reference range that is greater than that of Genova's healthy cohort.

↓↑ Indicates Genova's clinical condition cohort test results falling below and above the reference range that are greater than that of Genova's healthy cohort.

Cells with bolded arrows indicate Genova's clinical condition cohort had more test results falling above versus below ↓↑ or more below versus above ↓↑ the reference range compared to that of Genova's healthy cohort.

© Genova Diagnostics - Robert M. Davis, PhD, Lab Director - CLIA Lic. #11D025349 - Medicare Lic. #34-8475 - Georgia Lab Lic. Code #067-007  
New York Clinical Lab PFI #4578 - Florida Clinical Lab Lic. #800008124

- Features of this chart include:
  - A list of 16 biomarkers providing insight into digestive insufficiency, inflammation and gut microbiome metabolic activity
  - A column indicating any patient result that is out of reference range – either high (H) or low (L)
  - Clinical Associations for eight specific disease states

# Biomarker Clinical Association Chart

Patient: JANE DOE ID: Page 3

2200 GI Effects™ Comprehensive Profile – Stool

Biomarker	Patient Results Out of Reference Range	Genova Diagnostics Biomarker Clinical Associations*							
		IBS	IBD	Metabolic Syndrome	Chronic Fatigue	Auto-immune	Type 2 Diabetes	High Blood Pressure	Mood Disorders
Pancreatic Elastase		↓	↓	↓	↓	↓	↓	↓	↓
Products of Protein Breakdown (Total)	L						↓		
Fecal Fat (Total*)	H	↑		↑	↑	↑	↓	↑	↑
Triglycerides		↑			↑	↑	↓	↑	↑
Long Chain Fatty Acids	H	↑			↑	↑	↓	↑	↑
Cholesterol							↓	↑	
Phospholipids	H	↑	↑	↑	↑	↑	↑	↑	↑
Calprotectin			↑					↑	
Eosinophil Protein X (EPX)			↑						
Fecal sIgA		↑	↑	↑	↑	↑	↑	↑	↑
Short Chain Fatty Acids (SCFA) (Total)	L				↓	↓			
n-Butyrate Concentration	L			↓					
n-Butyrate %	H								
Acetate%	L				↓		↓		
Propionate %				↑			↑	↑	
Beta-glucuronidase						↓			↓

\*Information derived from GDx results comparing a healthy cohort to various clinical condition cohorts. The chart above showing a comparison of patient results to clinical conditions is meant for informational purposes only; it is not diagnostic, nor does it imply that the patient has a specific clinical diagnosis or condition.

The arrows indicate Genova's clinical condition cohort test results falling below or above the reference range that is greater than that of Genova's healthy cohort.

↑↓ Indicates Genova's clinical condition cohort test results falling below and above the reference range that are greater than that of Genova's healthy cohort.

Cells with bolded arrows indicate Genova's clinical condition cohort had more test results falling above versus below or more below versus above the reference range compared to that of Genova's healthy cohort.

© Genova Diagnostics - Robert M. Davis, PhD, Lab Director - CLIA Lic. #11D025349 - Medicare Lic. #34-0475 - Georgia Lab Lic. Code #067-007  
New York Clinical Lab PFI #4578 - Florida Clinical Lab Lic. #000008124

- Biomarkers found to have relevant associations with clinical conditions are displayed, along with the patient's results for these biomarkers
- These associations are:
  - Part of an ongoing research program
  - Will most likely evolve over time as data accrues and is refined
  - Reported for informational and clinical interest purposes only



# D

# I

# G

## GI Effects™ Comprehensive Profile - Stool

Methodology: GC/MS, Automated Chemistry, EIA

Results

QUINTILE DISTRIBUTION  
1st 2nd 3rd 4th 5th

Reference Range

### Digestion and Absorption

Pancreatic Elastase 1 †	>500		>200 mcg/g
Products of Protein Breakdown (Total*) (Valerate, Isobutylate, Isovalerate)	1.6 L		1.8-9.9 micromol/g
Fecal Fat (Total*)	59.4 H		3.2-38.6 mg/g
Triglycerides	1.2		0.3-2.8 mg/g
Long-Chain Fatty Acids	32.6 H		1.2-29.1 mg/g
Cholesterol	2.2		0.4-4.8 mg/g
Phospholipids	23.4 H		0.2-6.9 mg/g


### Inflammation and Immunology

Calprotectin †	<16		<=50 mcg/g
Eosinophil Protein X (EPX)†	0.4		<=4.6 mcg/g
Fecal secretory IgA	114		<=885 mcg/g

### Gastrointestinal Microbiome


<b>Metabolic</b>			
Short-Chain Fatty Acids (SCFA) (Total*) (Acetate, n-Butyrate, Propionate)	1.5 L		>=23.3 micromol/g
n-Butyrate Concentration	1.5 L		>=3.6 micromol/g
n-Butyrate %	100.0 H		11.8-33.3 %
Acetate %	<DL L		48.1-69.2 %
Propionate %	<DL		<=29.3 %
Beta-glucuronidase	2,393		368-6,266 U/g

The test results are grouped in three main areas: **D.**





# Pancreatic Elastase 1 (PE-1)



PE-1 is a proteolytic enzyme exclusively secreted by the pancreas

- Noninvasive biomarker of exocrine pancreatic function
- Reported according to Medical Decision Points

PE-1 ( $\mu$ /g)	Interpretation
200 to >500	Normal pancreatic exocrine function
100 to 199	Mild-to-Moderate insufficiency
<100	Severe insufficiency



# Pancreatic Exocrine Insufficiency (PEI)



World Journal of  
Gastroenterology

Online Submissions: <http://www.wjgnet.com/esps/bpgoffice@wjgnet.com>  
doi:10.3748/wjg.v19.i42.7258

World J Gastroenterol 2013 November 14; 19(42): 7258-7266  
ISSN 1007-9327 (print) ISSN 2219-2840 (online)  
© 2013 Baishideng Publishing Group Co., Limited. All rights reserved.

TOPIC HIGHLIGHT

Asbjørn Mohr Drewes, MD, PhD, DMSc, Professor, Series Editor

**Diagnosis and treatment of pancreatic exocrine insufficiency**

Björn Lindkvist

Björn Lindkvist, Institute of Medicine, Sahlgrenska Academy, University of Gothenburg, SE-413 45 Gothenburg, Sweden  
Björn Lindkvist, Division of Gastroenterology and Hepatology, Department of Internal Medicine, Sahlgrenska University Hospital, SE-413 45 Gothenburg, Sweden  
Author contributions: Lindkvist B designed and wrote the article.  
Correspondence to: Björn Lindkvist, MD, Associate Professor, Division of Gastroenterology and Hepatology, Department of Internal Medicine, Sahlgrenska University Hospital, SE-413 45 Gothenburg, Sweden. [bjorn.lindkvist@vgregion.se](mailto:bjorn.lindkvist@vgregion.se)  
Telephone: +46-31-3421000 Fax: +46-31-7412917  
Received: June 20, 2013 Revised: August 22, 2013  
Accepted: September 16, 2013  
Published online: November 14, 2013

spheres during meals. The dose should be in proportion to the fat content of the meal, usually 40-50000 lipase units per main meal, and half the dose is required for a snack. In cases that do not respond to initial treatment, the doses can be doubled, and proton inhibitors can be added to the treatment. This review focuses on current concepts of the diagnosis and treatment of pancreatic exocrine insufficiency.


© 2013 Baishideng Publishing Group Co., Limited. All rights reserved.

**Key words:** Chronic pancreatitis; Pancreatic exocrine insufficiency; Pancreatic enzyme replacement therapy

- Pancreatic exocrine insufficiency (PEI) is an important cause of maldigestion
- Symptoms in patients with PEI vary, depending on the degree and etiology of PEI
- Normal digestion requires adequate stimulation of pancreatic secretion, sufficient production of digestive enzymes by pancreatic acinar cells, and adequate mixing of pancreatic juice with ingested food
- The pancreatic juice plays a pivotal role in the digestion and absorption of nutrients



# Products of Protein Breakdown



## Products of Protein Breakdown (Isovalerate, Valerate and Isobutyrate)


- Derived from bacterial fermentation of protein in the distal colon
- Normal protein digestion and absorption is relatively complete in the stomach and small intestine

## Increased fecal Products of Protein Breakdown may be due to:

- Excessive delivery of protein to the colon, (ex: hypochlorhydria)
- Insufficient pancreatic proteases
- Excessive protein intake
- Increased colonic protein due to a GI bleed, excess mucus, bacterial overgrowth




# Fecal Fat (Total)



- The Fecal Fat test is a fecal fat extraction method that results in a quantitative value
- Elevated fecal fat may indicate maldigestion, malabsorption, or steatorrhea
  - Also query: transit time and dietary intake

The test results are grouped in three main areas: **D. I.**



# Calprotectin



- Neutrophil mediated inflammation in the gastrointestinal tract
  - Utilized for differentiating IBS from IBD
  - Quantify the degree of inflammation

Calprotectin (µg/g stool)	Interpretation	Follow up
≤50	Normal (no GI inflammation)	None
50 to 120	Borderline, suggestive of low grade inflammation	Re-evaluate at 4-6 weeks
>120	Abnormal	Determine source of inflammation and repeat test as clinically indicated
>250	Associated with high risk of clinical relapse	Adjust therapy accordingly

# Calprotectin


Original Article

## Does Fecal Calprotectin Predict Short-Term Relapse After Stopping Tnf $\alpha$ -Blocking Agents In Inflammatory Bowel Disease Patients In Deep Remission?

Pauliina Molander<sup>a,\*</sup>, Martti Färkkilä<sup>b,m</sup>, Ari Ristimäki<sup>c</sup>, Kimmo Salminen<sup>d</sup>, Helena Kemppainen<sup>d</sup>, Timo Blomster<sup>e</sup>, Ritva Koskela<sup>e</sup>, Airi Jussila<sup>f</sup>, Henna Rautiainen<sup>g</sup>, Markku Nissinen<sup>h</sup>, Johanna Haapamäki<sup>b</sup>, Perttu Arkkila<sup>b</sup>, Urpo Nieminen<sup>b</sup>, Juha Kuisma<sup>i</sup>, Jari Punkkinen<sup>j</sup>, Kaija-Leena Kolho<sup>k,m</sup>, Harri Mustonen<sup>l</sup>, Taina Sipponen<sup>b</sup>




*FC seems to increase and remain elevated before clinical or endoscopic relapse, suggesting that it can be used as a surrogate marker for predicting and identifying patients requiring close follow-up in clinical practice.*






# Eosinophil Protein X



- Eosinophil mediated inflammation in the gastrointestinal tract
  - Associated with inflammation and tissue damage
- Clinically, elevations in EPX indicate the presence of an IgE-mediated inflammatory process
  - Common associations are food allergies, parasitic infection, and IBD




# Fecal Secretory IgA



- Produced across mucosal surfaces, sIgA is the first line of defense in protecting the gastrointestinal epithelium from enteric toxins and pathogenic organisms
- Provides information on the competence of mucosal immunity and the potential risk for epithelial barrier dysfunction

The test results are grouped in three main areas: **D. I. G.**





# G/Gastrointestinal Microbiome

- Metabolic
  - Short Chain Fatty Acids
  - Beta glucuronidase
- Commensal Bacteria (PCR)
- Bacteriology & Mycology (culture and MALDI-TOF)
  - Sensitivities provided if applicable
- Parasitology (microscopy and EIA)
- Pathogens (culture and MALDI-TOF, add-on via EIA)




# Short Chain Fatty Acids

Metabolic					
Short-Chain Fatty Acids (SCFA) (Total*) (Acetate, n-Butyrate, Propionate)	1.5	L			$\geq 23.3$ micromol/g
n-Butyrate Concentration	1.5	L			$\geq 3.6$ micromol/g
n-Butyrate %	100.0	H			11.8-33.3 %
Acetate %	<DL	L			48.1-69.2 %
Propionate %	<DL				$\leq 29.3$ %

- Short Chain Fatty Acids (SCFA): Acetate, n-Butyrate and Propionate are produced by anaerobic bacterial fermentation of indigestible carbohydrate (fiber)
  - Important roles in maintaining colonic health and integrity
- Altered levels of fecal SCFA is an indicator of altered intestinal microbial composition



# Beta-glucuronidase




- Beta-glucuronidase is an enzyme induced by anaerobic bacteria
  - Many toxins, hormones, and drugs are excreted from the body after conjugation to a glucuronide molecule. Beta-glucuronidase can uncouple these conjugates, freeing these potential carcinogens in the bowel and promoting recirculation of these compounds through enterohepatic pathways
- A moderate level of Beta-glucuronidase activity is preferred as activity appears to be important for normal enterohepatic recirculation of endogenous compounds and vitamins
- Limited human studies show a relationship to colon cancer and hormone-related cancers
  - Calcium D-glucarate inhibits action of enzyme

# G/Gastrointestinal Microbiome

24 Commensal targets (PCR technology)


Commensal Balance and Relative Abundance graphics serve as broad markers of gut health



Commensal Bacteria (PCR)		Result CFU/g stool	QUINTILE DISTRIBUTION					Reference Range CFU/g stool
Bacteroidetes Phylum			1st	2nd	3rd	4th	5th	
<i>Bacteroides-Prevotella group</i>		1.4E9						3.4E6-1.5E9
<i>Bacteroides vulgatus</i>		4.8E9 H						<=2.2E9
<i>Barnesiella spp.</i>		4.3E8 H						<=1.6E8
<i>Odoribacter spp.</i>		2.3E8 H						<=8.0E7
<i>Prevotella spp.</i>		5.9E7 H						1.4E5-1.6E7
Firmicutes Phylum								
<i>Anaerotruncus colihominis</i>		4.3E7 H						<=3.2E7
<i>Butyrivibrio crossotus</i>		2.0E5						5.5E3-5.9E5
<i>Clostridium spp.</i>		2.6E10 H						1.7E8-1.5E10
<i>Coprococcus eutactus</i>		1.5E8 H						<=1.2E8
<i>Faecalibacterium prausnitzii</i>		1.2E10 H						5.8E7-4.7E9
<i>Lactobacillus spp.</i>		8.1E9 H						8.3E6-5.2E9
<i>Pseudoflavonifractor spp.</i>		4.3E8 H						4.2E5-1.3E8
<i>Roseburia spp.</i>		1.6E9						1.3E8-1.2E10
<i>Ruminococcus spp.</i>		5.2E7 L						9.5E7-1.6E9
<i>Veillonella spp.</i>		1.2E8 H						1.2E5-5.5E7
Actinobacteria Phylum								
<i>Bifidobacterium spp.</i>		3.2E9						<=6.4E9
<i>Bifidobacterium longum</i>		2.1E8						<=7.2E8
<i>Collinsella aerofaciens</i>		8.0E8						1.4E7-1.9E9
Proteobacteria Phylum								
<i>Desulfovibrio piger</i>		5.6E5						<=1.8E7
<i>Escherichia coli</i>		4.2E8 H						9.0E4-4.6E7
<i>Oxalobacter formigenes</i>		3.7E7 H						<=1.5E7
Euryarchaeota Phylum								
<i>Methanobrevibacter smithii</i>		2.0E8 H						<=8.6E7
Fusobacteria Phylum								
<i>Fusobacterium spp.</i>		1.8E4						<=2.4E5
Verrucomicrobia Phylum								
<i>Akkermansia muciniphila</i>		5.9E7						>=1.2E6




# Firmicutes/Bacteroidetes Ratio



- The F/B Ratio is a comparison of the 2 largest phyla
- F/B Ratio has been associated with disruption of metabolic homeostasis
  - Type 2 diabetes
  - Non-alcoholic fatty liver disease
  - Body Mass Index (BMI)



# Firmicutes/Bacteroidetes Ratio



**Fig. 1** The relative abundance of the major microbial phyla in different BMI categories (**a** BMI < 18.5, **b** BMI 18.5–24.9, **c** BMI 25–29.9 and **d** BMI ≥ 30)



## Gastrointestinal Microbiome

### Bacteriology (Culture)

*Lactobacillus spp.*

NG

*Escherichia coli*

4+ NP

*Bifidobacterium*

3+ NP

### Additional Bacteria

*Klebsiella pneumoniae*

4+ PP

*Bacillus species*

1+ NP


*gamma haemolytic Streptococcus*

1+ NP

### Mycology (Culture)

*Candida krusei*

1+ NP



### Microbiology Legend

NG



No Growth

NP



Non-Pathogen

PP



Potential Pathogen

P



Pathogen

# Bacteriology and Mycology

- Methodology is culture with MALDI-TOF MS identification
  - Vitex – MS (mass spectrometry) is an FDA-cleared platform for the rapid identification of bacteria and yeast from isolated culture colonies
  - Matrix Assisted Laser Desorption Ionization Time-Of-Flight
- MALDI-TOF examines patterns of proteins detected from intact bacteria and yeast isolated cultures
  - Technology ionizes chemical compounds (a combination of the culture sample to be examined and an assay matrix) to generate charged molecules, which provide a ‘tell-tale’ signature for identification





## Bacteria Sensitivity

### Prescriptive Agents

<i>Klebsiella pneumoniae</i>	R	I	S-DD	S	NI
Ampicillin	R				
Amox./Clavulanic Acid				S	
Cephalothin				S	
Ciprofloxacin				S	
Tetracycline				S	
Trimethoprim/Sulfa				S	

### Natural Agents

<i>Klebsiella pneumoniae</i>	LOW INHIBITION	HIGH INHIBITION
Berberine		
Oregano		
Plant Tannins		
Uva-Ursi		



## Mycology Sensitivity

### Azole Antifungals

<i>Candida krusei</i>	R	I	S-DD	S	NI
Fluconazole					32
Voriconazole				0.25	

### Non-absorbed Antifungals

<i>Candida krusei</i>	LOW INHIBITION	HIGH INHIBITION
Nystatin		

### Natural Agents

<i>Candida krusei</i>	LOW INHIBITION	HIGH INHIBITION
Berberine		
Caprylic Acid		
Garlic		
Undecylenic Acid		
Plant tannins		
Uva-Ursi		



# G/Gastrointestinal Microbiome – Parasitology

Parasitology

Microscopic Exam Results\*\*

No Ova or Parasites seen

Parasitology

Parasite Recovery: Literature suggests that >90% of enteric parasitic infections may be detected in a sample from a single stool collection. Increased sensitivity results from the collection of additional specimens on separate days.

Parasitology EIA Tests:

	In Range	Out of Range
<i>Cryptosporidium</i> ♦	Negative	
<i>Giardia lamblia</i> ♦	Negative	
<i>Entamoeba histolytica</i> ♦	Negative	



# Additional Results

Additional Results		
	Result	Expected Value
Fecal Occult Blood♦	Negative	Negative
Color††	Brown	
Consistency††	Formed/Normal	

Macroscopic Exam for Worms**
No larvae seen. Specimen contains mucus/mucosal material.




# Zonulin (add-on)

Stool Zonulin		
	Result	Reference Range
Zonulin, Stool	62.5	22.3-161.1 ng/mL

- Zonulin is a protein modulator of intestinal tight junctions (cellular structures that are critical for maintaining the integrity of the intestinal barrier)
- Release of zonulin has been shown to alter intestinal barrier integrity
- When tight junctions open, the intestinal barrier becomes permeable (intestinal permeability)
  - This is important for normal physiologic function
  - Prolonged release of zonulin is correlated with a dysfunctional intestinal barrier

# Zonulin (add-on)

- Prolonged compromise of intestinal barrier function (altered intestinal permeability) has been associated with a wide range of conditions encompassing gastrointestinal as well as systemic conditions
- Patients in whom zonulin testing may be most useful include:
  - Celiac Disease and other GI disorders (ex: non-celiac gluten sensitivity, IBS)
  - Cardiometabolic diseases including diabetes, obesity, NAFLD, insulin resistance, and associations with cardiometabolic risk and inflammatory markers






# Fecal Occult Blood

	Result	Expected Value
Fecal Occult Blood♦	Negative	Negative

- Identifies microscopic amounts of blood in the stool
- Hemosure diagnostic kit utilizes immunochemical technology (FIT)
- Higher sensitivity and specificity
  - No medication, supplement, or dietary restrictions
  - Specific to human hemoglobin (monoclonal and polyclonal antibodies)
  - One positive test warrants further investigation

# Case Example

- 52 y/o female with bloating and excessive belching and flatulence
- Constipation; uses a laxative three times per week to achieve a bowel movement
- Cerebellar ataxia
- Overweight
- Sedentary lifestyle
- WNL thyroid function
- Eats a standard American diet; sugar cravings & avoids fiber (aggravates symptoms)





# GI Effects® Stool Profile Overview & Considerations

- Maldigestion/absorption of fats
- Low dietary intake of meat; evaluate protein consumption
- Low SCFA; patient avoids fiber because it leads to an increase in gas and bloating
- PP Bacteria Klebsiella
- Elevated relative abundance compared to healthy cohort

**GI Effects profile can not identify SIBO, however there are markers that can be suggestive in the right patient population**



## Open

# Hydrogen and Methane-Based Breath Testing in Gastrointestinal Disorders: The North American Consensus

Ali Rezaie, MD, MSc, FRCP(C)<sup>1</sup>, Michelle Buresi, MD<sup>2</sup>, Anthony Lembo, MD<sup>3</sup>, Henry Lin, MD<sup>4</sup>, Richard McCallum, MD<sup>5</sup>, Satish Rao, MD<sup>6</sup>, Max Schmulson, MD<sup>7</sup>, Miguel Valdovinos, MD<sup>8</sup>, Salam Zakko, MD<sup>9</sup>, Mark Pimentel, MD, FRCP(C)<sup>1</sup> and on behalf of The North American Consensus group on hydrogen and methane-based breath testing

**OBJECTIVES:** Breath tests (BTs) are important for the diagnosis of carbohydrate maldigestion syndromes and small intestinal bacterial overgrowth (SIBO). However, standardization is lacking regarding indications for testing, test procedures, and interpretation. The purpose of this consensus was to develop recommendations for the interpretation of breath test results.

**METHODS:** Pre-meeting interpretation of breath test results was performed by a live meeting of experts in a live meeting, and a live meeting of experts in a live meeting, and a live meeting of experts in a live meeting.

**RESULTS:** Consensus was reached on the interpretation of breath test results. A rise in hydrogen of  $\geq 20$  ppm from baseline in hydrogen by 90 min should be considered a positive test for SIBO. Two peaks on breath test are not required for the diagnosis of SIBO. A level of  $\geq 10$  ppm be considered positive for methane on a breath test.

**CONCLUSIONS:** BT is a useful tool for the diagnosis of carbohydrate maldigestion syndromes and SIBO. The interpretation of breath test results should be based on the consensus recommendations.

**SUPPLEMENTARY MATERIAL** is linked to the online version of this article.

*Am J Gastroenterol* 2017; 112:775–784; doi:10.1038/ajg.2017.46; published online 21 March 2017


## INTRODUCTION

Breath test (BT) is performed to aid in the diagnosis of many common gastroenterological conditions including small intestinal bacterial overgrowth (SIBO) and irritable bowel syndrome

(IBS)-like symptoms, carbohydrate maldigestion and dysfunction or alterations in oro-cecal transit. Presently in clinical practice, BT is being performed with various substrates (e.g., glucose, lactulose, fructose, sorbitol, sucrose and inulin) using variable doses for a range

## Interpretation of breath testing results:

- “A rise of  $\geq 20$  ppm from baseline in hydrogen by 90 min should be considered a positive test for SIBO”
- “Two peaks on breath test are not required for the diagnosis of SIBO”
- “A level of  $\geq 10$  ppm be considered positive for methane on a breath test”



Evaluation for Methane (CH <sub>4</sub> )			
Peak methane level at any point			
	Result	Expected Value	
CH <sub>4</sub> Peak	>81	H	<10 ppm
A peak methane level ≥ 10 ppm at any point is indicative of a methane-positive result.			

Hydrogen (H <sub>2</sub> ), Methane (CH <sub>4</sub> ) and Carbon Dioxide								
	Baseline 0 min (S1)	20 min (S2)	40 min (S3)	60 min (S4)	90 min (S5)	120 min (S6)	150 min (S7)	180 min (S8)
H <sub>2</sub>	2	<2	<2	<2	<2	6	11	7
CH <sub>4</sub>	>81	>81	>81	>81	>81	50	78	64
H <sub>2</sub> + CH <sub>4</sub>	NR	NR	NR	NR	NR	56	89	71
CO <sub>2</sub> **	✓	✓	✓	✓	✓	✓	✓	✓
Actual Collection Times								
Actual Time	6:10	6:30	6:50	7:10	7:40	8:10	8:40	9:10
Actual Interval	0 min	20 min	40 min	60 min	90 min	120 min	150 min	180 min
**CO <sub>2</sub> is measured for quality assurance. ✓ indicates the CO <sub>2</sub> level is acceptable. ✗ indicates room air contamination exceeding acceptable limits.								




# Treatment

- Botanical Therapy for SIBO and PPs: garlic, goldenseal (berberine), and oregano x 4 weeks
- Iberogast, a prokinetic, was recommended
- Magnesium and exercise to assist with regularity
- Recommended increasing water intake to at least 1-2 L daily
- While fiber is important for constipation, it was avoided initially to address SIBO
- Encouraged increasing plant-based foods, healthy fats, and moderate amount of clean and lean meats
- Probiotics 4 hours away from botanical therapy

# Case Example #2

- 61 y/o female with IBS diagnosis
- History of bowel obstruction
- Bloating and alternating stools







Digestion and Absorption			
Pancreatic Elastase 1 †	452		>200 mcg/g
Products of Protein Breakdown (Total*)	4.4		1.2-9.9 micromol/g

## Gastrointestinal Microbiome

### Metabolic

Short-Chain Fatty Acids (SCFA) (Total\*)  
(Acetate, n-Butyrate, Propionate)

33.6



$\geq 23.3$  micromol/g

n-Butyrate Concentration

8.2



$\geq 3.6$  micromol/g

n-Butyrate %

24.4



11.8-33.3 %

Acetate %

56.6



48.1-69.2 %

Propionate %

18.9



$\leq 29.3$  %



Beta-glucuronidase

1,821



368-6,266 U/g

n-Butyrate %	24.4		11.8-33.3 %
Acetate %	56.6		48.1-69.2 %
Propionate %	18.9		$\leq 29.3$ %
Beta-glucuronidase	1,821		368-6,266 U/g



Gastrointestinal Microbiome							
Commensal Bacteria (PCR)	Result CFU/g stool	QUINTILE DISTRIBUTION					Reference Range CFU/g stool
		1st	2nd	3rd	4th	5th	
<b>Bacteroidetes Phylum</b>							
<i>Bacteroides-Prevotella</i> group	2.2E8						3.4E6-1.5E9
<i>Bacteroides vulgatus</i>	1.3E9						<=2.2E9
<i>Barnesiella</i> spp.	<DL						<=1.6E8
<i>Odoribacter</i> spp.	1.6E8 H						<=8.0E7
<i>Prevotella</i> spp.	7.1E6						1.4E5-1.6E7
<b>Firmicutes Phylum</b>							
<i>Anaerotruncus colihominis</i>	6.9E6						<=3.2E7
<i>Butyrivibrio crossotus</i>	<DL L						5.5E3-5.9E5
<i>Clostridium</i> spp.	1.7E9						1.7E8-1.5E10
<i>Coprococcus eutactus</i>	<DL						<=1.2E8
<i>Faecalibacterium prausnitzii</i>	1.3E9						5.8E7-4.7E9
<i>Lactobacillus</i> spp.	1.9E9						8.3E6-5.2E9
<i>Pseudoflavonifractor</i> spp.	7.1E7						4.2E5-1.3E8
<i>Roseburia</i> spp.	6.0E8						1.3E8-1.2E10
<i>Ruminococcus</i> spp.	<DL L						9.5E7-1.6E9
<i>Veillonella</i> spp.	4.1E6						1.2E5-5.5E7
<b>Actinobacteria Phylum</b>							
<i>Bifidobacterium</i> spp.	7.0E7						<=6.4E9
<i>Bifidobacterium longum</i>	9.0E6						<=7.2E8
<i>Collinsella aerofaciens</i>	6.7E7						1.4E7-1.9E9
<b>Proteobacteria Phylum</b>							
<i>Desulfovibrio piger</i>	<DL						<=1.8E7
<i>Escherichia coli</i>	2.3E7						9.0E4-4.6E7
<i>Oxalobacter formigenes</i>	4.9E6						<=1.5E7
<b>Euryarchaeota Phylum</b>							
<i>Methanobrevibacter smithii</i>	2.9E7						<=8.6E7
<b>Fusobacteria Phylum</b>							
<i>Fusobacterium</i> spp.	6.4E3						<=2.4E5
<b>Verrucomicrobia Phylum</b>							
<i>Akkermansia muciniphila</i>	<DL						>=1.2E6



## Gastrointestinal Microbiome

### Bacteriology

*Lactobacillus spp.*

*Escherichia coli*

*Bifidobacterium*

### Additional Bacteria

*alpha haemolytic Streptococcus*

*Mucoid Escherichia coli*

*Klebsiella pneumoniae*

*gamma haemolytic Streptococcus*

*Enterobacter cloacae*

### Mycology (Culture)

*Geotrichum species*

#### Bacteria Sensitivity

##### Prescriptive Agents

*Klebsiella pneumoniae*

R

I

S-DD

S

NI

Ampicillin

R

Amox./Clavulanic

Cephalothin

Ciprofloxacin

Tetracycline

Trimethoprim/Sulf

##### Natural Agents

*Klebsiella pneumoniae*

Berberine

Oregano

Plant Tannins

Uva-Ursi

##### Prescriptive Agents

*Enterobacter cloacae*

R

I

S-DD

S

NI

Ampicillin

R

Amox./Clavulanic Ac

Cephalothin

Ciprofloxacin

Tetracycline

Trimethoprim/Sulfa

##### Natural Agents

*Enterobacter cloacae*

Berberine

Oregano

Plant Tannins

Uva-Ursi

#### Mycology Sensitivity

##### Azole Antifungals

*Geotrichum species*

R

I

S-DD

S

NI

Fluconazole

32

Voriconazole

0.25

##### Non-absorbed Antifungals

*Geotrichum species*

LOW INHIBITION

HIGH INHIBITION

Nystatin

##### Natural Agents

*Geotrichum species*

LOW INHIBITION

HIGH INHIBITION

Berberine

Caprylic Acid

Garlic

Undecylenic Acid

Plant tannins

Uva-Ursi



## Parasitology

### Microscopic Exam Results\*\*

Blastocystis hominis: Many

### Parasites - Blastocystis spp. infection

*Blastocystis* spp.




Frequently Asked Questions (FAQs)

Biology

Resources for Health Professionals

Publications



*Blastocystis* is a common microscopic organism that inhabits the intestine and is found throughout the world. A full understanding of the biology of *Blastocystis* and its relationship to other organisms is not clear, but is an active area of research.

Image: B. spp. cyst-like forms in wet mounts under differential interference contrast (DIC) microscopy. Credit: [DPDx](#)



#### BLASTOCYSTIS SPP. FAQs

Most common questions answered...



#### PUBLICATIONS

Related articles and journal entries...



#### BIOLOGY

Stages of parasite development...



#### RESOURCES FOR HEALTH PROFESSIONALS

What you need to know...

### Parasitology EIA Tests:

In Range

Out of Range

*Cryptosporidium* ♦

Negative

*Giardia lamblia* ♦

Negative

*Entamoeba histolytica* ♦

Negative



# Treatment

- Parasite Treatment - Alinia BID x 7 days
  - Retest in 6 weeks
- Botanicals: Oil of oregano, horse chestnut, and berberine x 4 weeks
- Probiotics 4 hours away from therapy
- Decrease sugar, yeasted and molded foods
- Insoluble fiber and encourage food prebiotics (onion, garlic, artichoke, etc...)





**Lahnor Powell, ND, MPH**  
Moderator



**Michelle Maddux, ND**  
Presenter

Explore  
**WWW.GDX.NET**

*for more information and  
educational resources, including...*

LEARN **GDX** – Brief video modules

LIVE **GDX** – Previous webinar recordings

**GI University** – Focused learning modules

**Conferences** – Schedule of events we attend

**Test Menu** – Detailed test profile information

---

MY **GDX** – Order materials and get results

# Questions?





# Additional Questions?

**US Client Services: 800-522-4762**

**UK Client Services: 020-8336-7750**

**Please schedule a complimentary appointment with one of our Medical Education Specialists for questions related to:**

- Diagnostic profiles featured in this webinar
- How Genova's profiles might support patients in your clinical practice
- Review a profile that has already been completed on one of your patients

***We look forward to hearing from you!***



# Upcoming <sup>LIVE</sup> GDX Webinar Topics

January 24, 2018 – Susan Blum, MD, MPH Presents:

## *Arthritis: Taming the Flame by Treating Oxidative Stress*

Register for upcoming <sup>LIVE</sup> GDX Webinars online at [WWW.GDX.NET](http://WWW.GDX.NET)





**GI** *fx* **GI Effects**  
**Stool Profiles**



# The GI Effects® Advanced Interpretation

Digging Deeper

Michelle Maddux, ND

**GENOVA**  
DIAGNOSTICS®