



ESTRONEXSM PROFILE GUIDE



ANALYTES:

2-Hydroxyestrone (2-OHE1)
2-Hydroxyestradiol (2-OHE2)
4-Hydroxyestrone (4-OHE1)
16 α -Hydroxyestrone (16 α -OHE1)
2-Methoxyestrone (2-OMeE1)
4-Methoxyestrone (4-OMeE1)
Ratios
(2-OHE1+2-OHE2)/16 α -OHE1
2-OHE1/2-OMeE1

TESTING

The Estronex Profile is an easy-to-use and cost-effective way to assess estrogen metabolism. While estrogen-sensitive cancers are caused by multiple genetic and environmental factors, how effectively a woman metabolizes and eliminates estrogens gives information about her risk of developing certain cancers. Breast cancer is the second most common cause of cancer death in women.¹ As part of a breast cancer prevention program, women of all ages may benefit from the Estronex Profile. Women age 40 and older may consider taking the test regularly in addition to routine breast cancer screening methods. Breast cancer survivors may benefit from regular assessment of the Estronex Profile to prevent recurrence. The Estronex Profile can also be used to evaluate men at risk of developing breast and prostate cancer.^{2,3}

IMBALANCES IN ESTROGEN HAVE BEEN FOUND IN THE FOLLOWING CONDITIONS:⁴⁻¹⁵

- Breast cancer
- Colorectal cancer
- Prostate cancer
- Uterine cancers, including endometrial
- Ovarian cancer
- Cervical cancer
- Head and neck cancers
- Polycystic ovarian syndrome
- Uterine fibroids
- Endometriosis
- Preeclampsia
- Osteoporosis
- Menopausal symptoms
- Premenstrual syndrome

The Estronex Profile requires only a first morning urine collection. Research has shown that the 2/16 ratio is stable over the menstrual cycle but it is recommended that the test be completed at the same time in the menstrual cycle when comparing before- and after-treatment tests. After initiation of treatment, retesting is recommended within three to six months until the ratio is within normal limits, and then annually or biannually thereafter depending on degree of risk. Chemotherapeutics targeting estrogen or estrogen receptors may also change estrogen metabolism and monitoring the effects of such agents on estrogen levels may be useful.



BACKGROUND

Estrogen's basic biological function is as an anabolic hormone. There are three primary forms of estrogen:¹⁹ estrone (E1); estradiol (E2); and estriol (E3). All three forms have varying levels of activity and specificity for different tissues. E2 is the most active of the estrogens and is responsible for most of the actions attributed to the estrogens.

Because hormones exert great biological activity in small doses, it is very important that the body effectively eliminates them. Estrogens are cleared by phase I and phase II detoxification in the liver. In phase I, the estrogens are oxidized by the cytochrome P450 (CYP450) enzyme




FIGURE I. ESTROGEN METABOLISM

system (Figure 1). When CYP450 enzymes act on estradiol, the products are: 2-hydroxyestrone (2-OHE1), 4-hydroxyestrone (4-OHE1), and 16 α -hydroxyestrone (16 α -OHE1).¹⁶ Metabolism by CYP1A1 produces 2-OHE1, CYP1B1 produces 4-OHE1, and CYP2E family produces 16 α -OHE1.¹⁸ There are large variations in how individual women metabolize estrogen, which may be related to single nucleotide polymorphisms (SNPs) in these enzymes.



In women not on hormone-replacement therapy, 2-OHE1 is considered a “good” estrogen because it is associated with reduced cancer growth.¹⁷ 16 α -OHE1 has been shown to encourage tumor development and 4-OHE1 might cause DNA damage. Both 16 α -OHE1 and 4-OHE1 have been referred to as the “bad” estrogens.¹⁸⁻¹⁹ Poor ability to form 2-OHE1 leads to greater production of 16 α -OHE1.

In post-menopausal women on hormone-replacement therapy, a very high 2-OHE1 can result from high doses of estradiol. Thus, elevated levels of 2-OHE1 are associated with increased cancer risk.²² The prescribing doctor should be consulted about dosing when such very high 2-OHE1 are found.

INTERPRETATIONS

RATIOS

HYDROXYESTROGENS RATIO

Studies have shown that measuring the ratio of 2-hydroxylated estrogens (2-OHE1 + 2-OHE2) to 16 α -OHE1 (2/16 ratio) provides an important indication of risk for future development of breast cancer.²⁰ This ratio of “good to bad estrogen” can be determined from a single urine specimen. Studies have shown that individuals with breast cancer have lower 2/16 ratios. A low 2/16 ratio may indicate increased long-term risk for other estrogen-sensitive cancers including uterine, ovarian, cervical, prostate, and even head and neck cancers.^{18,21} Though a

woman’s biochemical individuality certainly affects how much of each form of estrogen is produced, in many women the primary cause of a low 2/16 ratio is related to dietary factors. Imbalances in the 2/16 ratio have been associated with prostate cancer,²³ colorectal cancer, and polycystic ovarian syndrome.²⁴

Low 2/16 Ratio

Urinary 2-OHE1 should be high enough relative to the urinary levels of 16 α -OHE1, to yield a 2/16 ratio that is greater than the low limit shown on the chart in the report.

Treatments To Increase 2/16 Ratio

- Decrease coffee consumption and smoking
- Lower an elevated BMI
- Supplement with indole-3-carbinol (I3C) or di-indolemethane (DIM)
- Increase cruciferous vegetables including: broccoli, cabbage, and Brussels sprouts
- Supplement with dried organic Brussels sprouts and kale
- Increase ground flax seed
- Increase intake of fish oils (omega-3 fatty acids)
- Increase fruit and vegetable intake
- Soy isoflavones (considered mildly estrogenic and may not be ideally suited for patients with hormone-sensitive cancers)
- Consider alternatives to oral contraceptive use (oral contraceptive users had significantly lower 2/16 ratios)

- Decrease inflammation and oxidative damage:
 - » Evaluate markers for oxidative stress, such as p-hydroxyphenyllactate, 8-hydroxy-2'-deoxyguanosine, and lipid peroxides
 - » Increase antioxidants in foods or supplements
- Ensure adequate estrogen detoxification and biotransformation:
 - » Sulfur-rich compounds such as N-acetyl cysteine, calcium D-glucarate and methyl donors such as S-adenosylmethionine (SAME), vitamin B12, folic acid, and vitamin B6
 - » Ensure effective GI elimination:
 - › Use 4R Protocol*
 - › Increase fiber, fruits, and vegetables
 - › Consume a low-allergen diet

* Refer to *Metametrix Handbook page 133*

High 2/16 Ratio

Most, but not all, epidemiological studies have found that women with a high 2/16 ratio are at reduced risk for breast cancer, however, excessive reduction of 16 α -OHE1 may reduce bone formation and increase oxidant stress.



Treatments to Decrease 2/16 Ratio

- Consider decreasing therapies that increase the 2/16 ratio (see above)
- Increase antioxidants
- Check methylation activity (homocysteine, methylmalonate, formiminoglutamate), accumulation of 2-OHE1 may demonstrate poor methylation activity
- Check bone mineral density and follow with the deoxypridinoline (DPD) test

METHYLATION RATIO

In phase II detoxification, the estrogen metabolites are conjugated with compounds that allow excretion by the body. Adequate methylation is very important to properly clear estrogens. High catechol-O-methyltransferase (COMT) activity has been associated with reduced risk for certain cancers. 2-Hydroxyestrone and 4-hydroxyestrone are methylated in Phase II to produce 2-methoxyestrone (2-OMeE1) and 4-methoxyestrone (4-OMeE1). 2-OMeE1 has anticancer effects and studies show that it is lower in breast cancer patients than controls. It has been referred to as a “good” estrogen. 4-OMeE1 is a non-carcinogenic metabolite of 4-OHE1.

High 2-OHE1/2-OMeE1 Ratio

- Imbalanced estrogen metabolism
- Low activity COMT
- Evaluate methylation activity:
 - » COMT genetic testing
 - » Homocysteine
 - » B₁₂
 - » Folic Acid
 - » SAME
 - » B₆
 - » Reduce stress: COMT is involved in the metabolism of epinephrine, reducing availability for estrogen metabolism

ANALYTES

While the ratios are important, the clinician should note the values of the individual estrogens. Ideally, the patient's levels of individual estrogen metabolites will be within normal limits for the population group (premenopausal, post-menopausal without hormone therapy, post-menopausal with hormone therapy, or males).



2-HYDROXYESTRONE (2-OHE1)

High 2-OHE1

- High levels of the “good estrogen” are typically considered beneficial; however, when coupled with a low 2-OMeE1, it may indicate poor methylation activity (see high 2-OHE1/2-OMeE1).
- If both 2-OHE1 and 16 α -OHE1 are high, then total estrogen load may be high. Consider methods to better clear the estrogens from the body or HRT dosing may need revision if only 2-OHE1 is very high.

Low 2-OHE1

- Low levels of the “good estrogen” (see treatments to increase 2/16 ratio on page 3)
- If 2-OHE1 is low and 16 α -OHE1 is low, then total estrogen production may be low. Address underlying causes of low hormone production, as needed.

16 ALPHA-HYDROXYESTRONE (16 α -OHE1)

High 16 α -OHE1

- High levels of the “bad estrogen” (see treatments to increase 2/16 ratio)
- If both 2-OHE1 and 16 α -OHE1 are high, then total estrogen load may be high. Consider methods to better biotransform and excrete the estrogens from the body or decrease HRT.

Low 16 α -OHE1

- Low levels of the “bad estrogen” are typically considered beneficial.
- If 2-OHE1 is low and 16 α -OHE1 is low, then total estrogen production may be low. Address underlying causes of low hormone production.

2-METHOXYESTRONE (2-OMeE1)

2-OMeE1 is produced from 2-OHE1 through the COMT enzyme. Anti-cancerogenic effects have been ascribed to 2-OHE1 and particularly 2-OMeE1.²³ 2-OMeE1 has shown antiproliferative effects in both hormone-dependent and hormone-independent breast cancer cells.^{25,26} These studies have shown that urinary 2-OMeE1 levels were lower in breast cancer patients than controls.

High 2-OMeE1

- Generally no treatment recommended
- Ensure adequate phase II detoxification

Low 2-OMeE1

Identify cause:

- If 2-OHE1 and 4-OHE1 are high, then low 2-OMeE1 is likely the result of poor methylation. Follow methylation treatment recommendations.
- If 2/16 ratio is also low, then low 2-OMeE1 may be indicating a CYP imbalance. Follow recommendations for improving a low 2/16 ratio.
- If total estrogen metabolite production is low, then consider direct assessment of estrogen levels.
- Ensure adequate COMT function via methyl donors such as SAmE, B₁₂, folic acid, and B₆.
- Poor cancer protection
- Associated with preeclampsia

4-HYDROXYESTRONE (4-OHE1)

4-OHE1 is referred to as the “bad” estrogen, along with 16 α -OHE1. They are primarily produced by CYP1B1 and CYP34A, respectively, enzymes localized in tissues, including breast and prostate as well as liver. Some have suggested that increased expression of CYP1B1 and 4-hydroxylation of estradiol are biomarkers of tumorigenesis.^{27,28} 4-OHE1 may be further metabolized to DNA-damaging quinone estrogens.²⁹ It may also have estrogenic effects.³⁰ Human breast cancer tissue produces much higher levels of 4-OHE1 than 2-OHE1, while normal breast tissue produces approximately equal amounts of the two metabolites. Women taking hormone therapy with a polymorphism in CYP1B1 had twice the risk of developing breast cancer compared to other HRT users.³¹

High 4-OHE1

- High levels of the “bad estrogen”
- Improve methylation by adding cofactors (B₁₂, folate) or methyl donors (betaine, dimethyl glycine)
- Consider genetic testing for COMT and CYP1B1 activity, particularly if positive family history
- Reduce stress: COMT is involved in the metabolism of epinephrine, reducing availability for estrogen metabolism
- Increase inhibitors of CYP1B1
 - » Grapefruit
 - » Ginseng
- See treatment for 2/16 ratio
- Avoidance of CYP1B1 inducers
 - » Polycyclic aromatic hydrocarbons
- Evaluate methylation activity:
 - » Serum Homocysteine
 - » Serum B₁₂ or methylation
 - » Urinary FIGLU
 - » Urinary xanthurenate
- Other potential protective factors:
 - » Glutathione (reduction of estrogen quinones)
 - » Resveratrol (prevents estrogen quinone formation)
 - » Selenium, zinc, magnesium

Low 4-OHE1

- Low levels of the “bad estrogen”
- Generally no treatment recommended

4-METHOXYESTRONE (4-OMeE1)

4-OMeE1 is the COMT-detoxified, non-carcinogenic product of 4-OHE1.²⁹ Low levels of 4-OMeE1 were found in the urine of both high-risk and breast cancer patients.³² It may be protective to have higher levels of the deactivated product, 4-OMeE1.

High 4-OMeE1

- Generally no treatment recommended

Low 4-OMeE1

- Improve methylation, see treatments for high 4-OHE1

**OTHER TESTING**

There are other laboratory tests that reveal important information about a patient’s estrogen metabolism. Higher levels of inflammation can impact estrogen receptor reactivity.³³ Oxidative stress should be assessed to assure adequate hormone function. Markers for oxidative stress are p-hydroxyphenyllactate, 8-hydroxy-2'-deoxyguanosine, high sensitivity C-reactive protein (hs-CRP), and lipid peroxides. Tests to evaluate levels of anti-inflammatory fatty acids can also be useful. Hormones excreted in bile can be deconjugated by intestinal bacteria and reabsorbed in cases of gastrointestinal dysbiosis. If dysbiosis is suspected, consider a GI EffectsSM Profile, OrganixSM Dysbiosis Profile, and an AllergixSM Food Antibody Profile. Evaluating xenoestrogen levels (Phthalates & Parabens Profile) or single nucleotide polymorphisms (SNPs) in detoxification enzymes could be helpful when treatments to improve estrogen clearance are ineffective. Finally, a full urinary hormone profile may be indicated to evaluate a patient’s levels of other hormones.

REFERENCES

1. C.D.C. Cancer Prevention and Control. <http://www.cdc.gov/cancer/dcpc/data/women.htm>. Accessed July 16, 2009.
2. Paola Muti, Kim Westerlind, Tiejian Wu, et al. Urinary estrogen metabolites and prostate cancer: a case-control study in the United States. *Cancer Causes and Control*. Dec 2002; 13(10): 1573-7225.
3. Barba M, Yang L, Schünemann HJ, et al. Urinary estrogen metabolites and prostate cancer: a case-control study and meta-analysis. *J Exp Clin Cancer Res*. Oct 8 2009; 28:135.
4. Yoo HJ, Sepkovic DW, Bradlow HL, Yu GP, Sirilian HV, Schantz SP. Estrogen metabolism as a risk factor for head and neck cancer. *Otolaryngol Head Neck Surg*. Mar 2001;124(3):241-247.
5. Compston J. Clinical and therapeutic aspects of osteoporosis. *Eur J Radiol*. Aug 4 2009.
6. Lopez LM, Kaptein AA, Helmerhorst FM. Oral contraceptives containing drospirenone for premenstrual syndrome. *Cochrane Database Syst Rev*. 2009(2):CD006586.
7. Leyendecker G, Wildt L, Mall G. The pathophysiology of endometriosis and adenomyosis: tissue injury and repair. *Arch Gynecol Obstet*. Jul 31 2009.
8. Okolo S. Incidence, aetiology and epidemiology of uterine fibroids. *Best Pract Res Clin Obstet Gynaecol*. Aug 2008;22(4):571-588.
9. Veerapaneni P, Kirma N, Nair HB, Hammes LS, Hall KL, Tekmal RR. Elevated aromatase expression correlates with cervical carcinoma progression. *Gynecol Oncol*. Sep 2009;114(3):496-500.
10. Berstein L, Zimarina T, Imyanitiv E, et al. Hormonal imbalance in two types of endometrial cancer and genetic polymorphism of steroidogenic enzymes. *Maturitas*. Jul 20 2006;54(4):352-355.
11. Sherman ME. Theories of endometrial carcinogenesis: a multidisciplinary approach. *Mod Pathol*. Mar 2000;13(3):295-308.
12. Hopkinson ZE, Sattar N, Fleming R, Greer IA. Polycystic ovarian syndrome: the metabolic syndrome comes to gynaecology. *Bmj*. Aug 1 1998;317(7154):329-332.
13. Stenchever MA. *Comprehensive Gynecology*. 4th ed. St. Louis, MO: Mosby; 2001.
14. Ferreira AM, Westers H, Albergaria A, Seruca R, Hofstra RM. Estrogens, MSI and Lynch syndrome-associated tumors. *Biochim Biophys Acta*. Jun 25 2009.
15. Iqbal J, Zaidi M. Understanding estrogen action during menopause. *Endocrinology*. Aug 2009;150(8):3443-3445.
16. Cribb AE, Knight MJ, Dryer D, et al. Role of polymorphic human cytochrome P450 enzymes in estrone oxidation. *Cancer Epidemiol Biomarkers Prev*. Mar 2006;15(3):551-558.
17. Lord RS, Bongiovanni B, Bralley JA. Estrogen metabolism and the diet-cancer connection: rationale for assessing the ratio of urinary hydroxylated estrogen metabolites. *Altern Med Rev*. Apr 2002;7(2):112-129.
18. Muti P, Bradlow HL, Micheli A, et al. Estrogen metabolism and risk of breast cancer: a prospective study of the 2/16alpha-hydroxyestrone ratio in premenopausal and postmenopausal women. *Epidemiology*. 2000;11(6):635-640.
19. Lemon HM, Heidel JW, Rodriguez-Sierra JF. Increased catechol estrogen metabolism as a risk factor for nonfamilial breast cancer. *Cancer*. Jan 15 1992;69(2):457-465.
20. Kabat GC, O'Leary ES, Gammon MD, et al. Estrogen metabolism and breast cancer. *Epidemiology*. Jan 2006;17(1):80-88.
21. Yadav BS, Sharma SC, Patel FD, Ghoshal S, Kapoor R, Kumar R. Nonbreast second malignancies after treatment of primary breast cancer. *Int J Radiat Oncol Biol Phys*. Apr 1 2009;73(5):1489-1492.
22. Wellejus A, Olsen A, Tjonneland A, Thomsen BL, Overad K, Loft S. Urinary Hydroxyestrogens and Breast Cancer Risk among Postmenopausal Women: A Prospective Study. Institute of Public Health, University of Copenhagen, DE, and Dept. of Clinical Epidemiology, Aalborg Hosp., Aarhus Univ. Hosp. Aalborg, DE.
23. Hamilton-Reeves JM, Rebello SA, Thomas W, Slaton JW, Kurzer MS. Soy protein isolate increases urinary estrogens and the ratio of 2/16 alpha-hydroxyestrone in men at high risk of prostate cancer. *J Nutr*. Oct 2007;137(10):2258-2263.
24. Salih S, Xu X, Veenstra TD, et al. Lower levels of urinary 2-hydroxyestrogens in polycystic ovary syndrome. *J Clin Endocrinol Metab*. Aug 2007;92(8):3285-3291.
25. Singh A, Purohit A, Hejaz HA, Potter BV, Reed MJ. Inhibition of deoxyglucose uptake in MCF-7 breast cancer cells by 2-methoxyestrone and 2-methoxyestrone-3-O-sulfamate. *Mol Cell Endocrinol*. Feb 25 2000;160(1-2):61-66.
26. Thibodeau PA, Kachadourian R, Lemay R, Bisson M, Day BJ, Paquette B. In vitro pro- and antioxidant properties of estrogens. *J Steroid Biochem Mol Biol*. Jul 2002;81(3):227-236.
27. Spink BC, Fasco MJ, Gierthy JF, Spink DC. 12-O-tetradecanoylphorbol-13-acetate upregulates the Ah receptor and differentially alters CYP1B1 and CYP1A1 expression in MCF-7 breast cancer cells. *J Cell Biochem*. Sep 1 1998;70(3):289-296.
28. Hayes CL, Spink DC, Spink BC, Cao JQ, Walker NJ, Sutter TR. 17 beta-estradiol hydroxylation catalyzed by human cytochrome P450 1B1. *Proc Natl Acad Sci USA*. Sep 3 1996;93(18):9776-9781.
29. Yang L, Gaikwad NW, Meza J, et al. Novel biomarkers for risk of prostate cancer: results from a case-control study. *Prostate*. Jan 1 2009;69(1):41-48.
30. Castagnetta LA, Granata OM, Traina A, et al. Tissue content of hydroxyestrogens in relation to survival of breast cancer patients. *Clin Cancer Res*. Oct 2002;8(10):3146-3155.
31. Rylander-Rudqvist T, Wedren S, Granath F, et al. Cytochrome P450 1B1 gene polymorphisms and postmenopausal breast cancer risk. *Carcinogenesis*. Sep 2003;24(9):1533-1539.
32. Gaikwad NW, Yang L, Muti P, et al. The molecular etiology of breast cancer: evidence from biomarkers of risk. *Int J Cancer*. May 1 2008;122(9):1949-1957.
33. Tamir S, Izrael S, Vaya J. The effect of oxidative stress on ERalpha and ERbeta expression. *J Steroid Biochem Mol Biol*. Aug 2002;81(4-5):327-332.



3425 Corporate Way Duluth, GA 30096 800-221-4640 Fax 770-441-2237 www.metamatrix.com