

Step 3:

Returning The Samples To The Laboratory

Samples must be returned to Genova Diagnostics by overnight delivery to arrive next day.

For UK mainland patients only: Genova Diagnostics recommends using the enclosed courier collection service, or alternative e.g.: Royal Mail, Special Delivery. The mailing envelope is not pre-paid.


For patients outside of the UK: please arrange an international courier such as FedEx, DHL, UPS.

Please note: the laboratory is closed at weekends. Samples should not be in transit over the weekend; please arrange for the samples to be sent overnight and to arrive no later than Friday.

- Plan to return the samples **Monday – Thursday by overnight delivery only.**
- Make sure the collection tubes are tightly closed and identified with completed labels. Failure to complete the labels correctly will result in the samples being rejected.
- Seal the samples in the biohazard bags with the absorbent pad.
- Place the Day 1, Day 2 and Day 3 biohazard bags inside the mailing envelope with your completed and signed requisition form.
- Seal the envelope.
- Once posting/couriering of samples has been arranged, keep any tracking numbers issued by the courier / Royal Mail for future reference.

Parasitology / Macroscopic Worm Analysis

Patient Instructions



Check Your Kit

- A - 1 Requisition Form (to be completed and signed)
- B - 1 mailing envelope

Day One Bag

- C - 1 Green-top SAF tube
- D - 1 Collection tray
- E - 1 Glove
- F - 1 Absorbent pad
- G - 1 Label

Day Two Bag

- H - 1 Green-top SAF tube
- I - 1 Collection tray
- J - 1 Glove
- K - 1 Absorbent pad
- L - 1 Label

Day Three Bag

- M - 1 Green-top SAF tube
- N - 1 Collection tray
- O - 1 Glove
- P - 1 Absorbent Pad
- Q - 1 Label

- If any of the items are missing or expired, call the laboratory on 020 8336 7750.
- The appearance of the above kit components are for illustration use only and actual contents may vary slightly from that shown.

Step 1:

Important things to know and consider

- If any liquid is spilled from the vials, stop the test immediately. Any liquid lost from a vial **invalidates** the test results. **Call the laboratory on +44 (0) 20 8336 7750** to obtain new vial(s).
- **Extreme care** should be used when opening and filling vials. Do not spill or discard any tube fluid. Avoid contact of the eyes or skin with the fluid in the tube. For contact with eyes, wash thoroughly for 15 minutes. For skin contact, wash thoroughly with soap and water. Do not inhale or ingest liquid. For accidental ingestion, contact your healthcare provider at once.
- The following are **interfering substances** and should be discontinued **unless directed otherwise by your practitioner**: bentonite clay, betaine HCL or digestive enzymes, castor oil, mineral oil, rectal suppositories. If you are taking a medication that isn't listed as an interfering substance, ask your healthcare professional before stopping any prescribed substance.
- We do not recommend performing this test when stool samples are excessively dilute or watery.
- Wait at least **4 weeks** from colonoscopy or barium enema before beginning the test.
- **2 weeks before test**: If taking antibiotics or antifungals it is recommended that you **wait at least 14 days** after your last dose before beginning the test (unless instructed otherwise by your practitioner).
- **Eat your usual diet** during the 2-3 days prior to stool collection.

Schedule & Prepare for Stool Collection

- **The test is collected from 3 bowel movements. Prepare your specimen using these collection tubes on these days:**
 - Day One:** 1 Green-top SAF tube
 - Day Two:** 1 Green-top SAF tube
 - Day Three:** 1 Green-top SAF Tube

If unable to collect on consecutive days, resume the test on the next day you have a bowel movement.

- **Samples must be received within 4 days of Day 3 collection.** Plan to collect the samples and complete the test no later than Thursday, ensuring samples **arrive at the laboratory by Friday**. If you complete the test over the weekend, arrange to return the samples on Monday for Tuesday.
- **Storing your samples:** Send in your samples as soon as you complete the test. Store samples at room temperature or refrigerate. **Do NOT freeze.**
- **Complete the Requisition Form** with all patient and payment information. Be sure it is signed by the Patient/Responsible party in the box labelled final sample date / time.

Step 2:

Stool Collection

Not following these instructions may affect the test results.



- 1** Pass urine into the toilet if necessary. Pass stool into the collection tray. Do not let urine or water from toilet touch the stool specimen.



- 2** Open the green-top SAF collection tube. Using the spoon in the cap of the tube, pick up several portions of stool from contrasting areas in the tray. Be sure to collect any blood or pus and samples from areas of different colour or texture. Add sufficient stool sample for stool and tube liquid to reach the red "fill" line, as shown on the label. **Do not overfill.**



- 3** Mix and mash the specimen with the spoon until the sample is thoroughly mixed and is as smooth as possible.



- 4** Recap the tube and check to ensure it is securely capped. Then, shake well until mixed thoroughly.



- 5** Write your full name (first and last), date of birth, time and date of collection on the label enclosed in the bag and affix it to the tube.



- 6** Be certain the caps are tight before returning. Place your filled tube in the corresponding day's biohazard bag with the absorbent pad. Seal the bag securely. **Do NOT freeze** the sample.

- 7** On Day Two and Day Three, repeat steps 1 to 6

Consult your healthcare provider if you have any questions at any time during this test.